

# **SITEK RESEARCH LABORATORIES**

15235 Shady Grove Road, Suite 303, Rockville, Maryland 20850 • 301/926-4900 FAX 301/926-8891

## **FINAL REPORT**

### **Study Title**

Evaluation of a Test Article in the *Salmonella typhimurium*/*Escherichia coli*  
Plate Incorporation Mutation Assay in the Presence and Absence  
of Induced Rat Liver S-9

### **Test Article**

Ethylenediamine dinitrate (EDDN)

### **Authors**

Jian Song, Ph.D.  
Paul E. Kirby, Ph.D.

### **Performing Laboratory**

SITEK Research Laboratories  
15235 Shady Grove Road, Suite 303  
Rockville, Maryland 20850

### **Laboratory Project I.D.**

SITEK Study No.: 1003-2140

### **Study Initiation Date**

August 5, 2009

### **Study Completion Date**

February 25, 2010

### **Sponsor**

USA RDECOM, AMSRD-MSF  
Environmental Acquisition & Logistics Sustaining Program  
Aberdeen Proving Ground, MD 21010

### **Sponsor's Study Coordinator**

Gunda Reddy, Ph.D., DABT

Report Documentation Page				Form Approved OMB No. 0704-0188		
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.						
1. REPORT DATE <b>25 FEB 2009</b>		2. REPORT TYPE		3. DATES COVERED		
4. TITLE AND SUBTITLE <b>Evaluation of a Test Article in the Salmonella typhimurium/Escherichia coli Plate Incorporation Mutation Assay in the Presence and Absence of Induced Rat Liver S-9. Test Article - Diethylene diamine dinitrate (EDDN)</b>				5a. CONTRACT NUMBER <b>W91ZLK-09-P-0958</b>		
				5b. GRANT NUMBER		
				5c. PROGRAM ELEMENT NUMBER		
6. AUTHOR(S) <b>Jian Song; Paul Kirby</b>				5d. PROJECT NUMBER		
				5e. TASK NUMBER		
				5f. WORK UNIT NUMBER		
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) <b>SITEK Research Laboratories,15235 Shady Grove Road, Suite 303,Rockville,MD,20850</b>				8. PERFORMING ORGANIZATION REPORT NUMBER <b>1003-2140</b>		
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)		
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)		
12. DISTRIBUTION/AVAILABILITY STATEMENT <b>Approved for public release; distribution unlimited.</b>						
13. SUPPLEMENTARY NOTES						
14. ABSTRACT <b>Ethylenediamine dinitrate (EDDN) (99.5% pure) was tested for mutagenic potential with Salmonella typhimurium strains, TA 98, TA 100, TA 1535, TA 1537 and Escherichia coli strain WP2 uvrA by plate incorporation method according to OECD TG 471 in compliance with Good Laboratory Practice, EDDN was tested at doses 100, 500, 1000, 3000 ug/plate and 5000 ug/plate. The results showed that EDDN was not mutagenic in Salmonella typhimurium and Escherichia coli both with and without activation in tester strains (TA 98, TA 1535, and TA 1537) but it showed positive increase in TA 100 tester strain without activation (3.1-fold) at 5000 ug/plate. This was also verified with confirmatory mutation assay at concentrations up to 5000 ug/plate with activation and also found positive (2.5 fold) in TA 100. No positive increases were observed with other tester strains (TA 98, TA 1535, TA 1537) with or without activation. Based on the positive response was only in one strain TA 100 observed at the highest dose under conditions of test and according to criteria set for evaluating test results, the EDDN was weak positive in this system both with and without activation.</b>						
15. SUBJECT TERMS						
16. SECURITY CLASSIFICATION OF:				17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES <b>73</b>	19a. NAME OF RESPONSIBLE PERSON
a. REPORT <b>unclassified</b>	b. ABSTRACT <b>unclassified</b>	c. THIS PAGE <b>unclassified</b>				

## STUDY DIRECTOR'S COMPLIANCE STATEMENT

Study No.: 1003-2140

Sponsor's Test Article I.D.: Ethylenediamine dinitrate (EDDN)

The protocol for this study was designed to meet or exceed the US EPA, OECD, and ICH Guidelines specified in the following documents:

United States Environmental Protection Agency, Title 40 Code of Federal Regulations, Part 798, Health Effects Testing Guidelines, Subpart F, Sec. 798.5265, the *Salmonella typhimurium* reverse mutation assay. Revised July 1, 2002.

OECD Guideline for the Testing of Chemicals, No. 471. Bacterial Reverse Mutation Test. Revised July 21, 1997.

International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use. ICH Harmonized Tripartite Guideline S2A. Guidance on Specific Aspects of Regulatory Genotoxicity Tests for Pharmaceuticals. Federal Register 61 (80):18198-18202, 1996.

The study described in this report was conducted in compliance with the following listed Good Laboratory Practice standards with the exception that the dosing solutions analysis was not conducted:

United States Environmental Protection Agency, Title 40 Code of Federal Regulations Parts 160 and 792, Revised July 1, 2002.

United States Food and Drug Administration, Title 21 Code of Federal Regulations Part 58, Revised April 1, 2003.

Japanese Ministry of Agriculture, Forestry and Fisheries, 11 NohSan, Notification No. 6283, October 1, 1999.

Japanese Ministry of Health and Welfare, Ordinance No. 21, April 1, 1997.

Japanese Ministry of International Trade and Industry, Notification No. 85, Basic Industries Bu, March 31, 1984.

Organization for Economic Cooperation and Development, The OECD Principles of Good Laboratory Practice, Environment Monograph No. 45 [ENV/MC/CHEM(98)17], Paris 1998.

Signature Paul E. Kirby  
Paul E. Kirby, Ph.D.<sup>1</sup>  
Study Director

2-25-10  
Date

---

<sup>1</sup> Dr. Jian Song was the Study Director for the in-life phase of this study and was the author of the draft report. He was not in the employ of SITEK Research Laboratories when this final report was prepared, therefore, Dr. Kirby has replaced him as Study Director.

**QUALITY ASSURANCE UNIT'S STATEMENT**

Study No.: 1003-2140

Sponsor's Test Article I.D.: Ethylenediamine dinitrate (EDDN)


The performance of this study was audited for adherence to the Good Laboratory Practice regulations for nonclinical laboratory studies by the Quality Assurance Unit of SITEK Research Laboratories. In this context, the facilities, equipment, personnel, methods, practices, controls, original data and reports have been inspected as per SITEK's Quality Assurance Unit's Standard Operating Procedures. The information contained within this report accurately reflects the raw data generated from this study.

Protocol Review Date: 08-05-08

The following phases were inspected for this study:

<u>Inspection Date</u>	<u>Phases Inspected</u>	<u>Date Findings Reported to Study Director</u>	<u>Date Findings Reported to Management</u>
<u>08-27-09</u>	<u>Plating of Tester Strain Genotypes</u>	<u>08-31-09</u>	<u>09-01-09</u>
<u>09-18-09</u>	<u>Workbook Audit</u>	<u>09-22-09</u>	<u>09-23-09</u>
<u>09-23-09</u>	<u>Draft Report Audit</u>	<u>09-24-09</u>	<u>09-28-09</u>
<u>02-25-10</u>	<u>Final Report Audit</u>	<u>02-25-10</u>	<u>02-25-10</u>

Signature \_\_\_\_\_

  
 Ray Kirby, B.S.  
 Quality Assurance


2/25/10  
 Date


**STUDY DIRECTOR'S SIGNATURE PAGE**

This study was performed under the supervision of Jian Song, Ph.D., Study Director for *Salmonella typhimurium* and *Escherichia coli* Gene Mutation Assays, at SITEK Research Laboratories, 15235 Shady Grove Road, Suite 303, Rockville, Maryland 20850.

The Draft Report for this study was written by Dr. Song and released on September 28, 2009. The Final Report was prepared by Dr. Paul E. Kirby<sup>2</sup> and released on February 25, 2010.

**Signature**

  
\_\_\_\_\_  
Paul E. Kirby, Ph.D.  
Study Director

  
\_\_\_\_\_  
Date

---

<sup>2</sup> Dr. Song was no longer in the employ of SITEK Research Laboratories when the final report was prepared, therefore, Dr. Kirby has replaced him as Study Director.

## ABSTRACT

The test article, Ethylenediamine dinitrate (EDDN) (Lot NO.: ABY08L031S010, 99.5% pure), was tested for its potential to cause mutations at the histidine operon of *Salmonella typhimurium* strains TA98, TA100, TA1535 and TA1537, and at the tryptophan operon of *Escherichia coli* strain WP2 uvrA. The assay was conducted using the plate incorporation method of treatment.

The test was conducted according to the method of Ames et al. in the presence and absence of metabolic activation using the S-9 fraction prepared from livers of Aroclor 1254-induced rats. The test article was tested for toxicity to strains TA100 and WP2 uvrA in a Range Finding Test at concentrations ranging from 5.0-5000 µg/plate. The test article was dissolved and subsequently diluted in deionized, distilled water. Water was used as a solvent control. The tester strains were exposed to the test article in the absence of exogenous activation and in the presence of Aroclor 1254-induced rat liver S-9 plus cofactors. Toxicity was evaluated based on: 1) reversion frequency, 2) viability, and 3) integrity of the background lawn.

The results of the Range Finding Test for TA100 and WP2 uvrA indicated that the test article, EDDN was found not toxic to both TA100 and WP2 uvrA at all dose levels both with and without activation. The background lawn was normal.

Since no toxicity was observed in the Range Finding Test, the Definitive Ames Assay was conducted with the maximum dose level recommended in the regulatory guidelines as the top dose (5000 µg/plate) for *Salmonella typhimurium* and *Escherichia coli* both with and without activation. The other doses were 100, 500, 1000 and 3000. The results indicate that under the condition of this study, EDDN, did cause positive increases in the mean number of revertants per plates with the tester strain TA100 without activation (3.1-fold) at 5000 µg/plate. No positive increases were observed with other tester strain/activation combinations. The background lawns were normal. The solvent controls and positive controls fulfilled the requirements of a valid test.

To verify the result from the Definitive Mutation Assay a Confirmatory Mutation Assay was performed using the plate incorporation method at concentrations of 3000, 3500, 4000, 4500 and 5000 µg/plate for *Salmonella typhimurium* and *Escherichia coli* both with and without activation. The results confirm the finding of the Definitive Mutation Assay. EDDN, did cause positive increases in the mean number of revertants per plates with the tester strain TA100 without activation (2.6-fold) at 5000 µg/plate, also with activation was positive (2.5-fold) at 5000 µg/plate in the Confirmatory Mutation Assay. No positive increases were observed with other tester strain/activation combinations. The background lawns were normal. The solvent controls and positive controls fulfilled the requirements of a valid test.

The results of the Mutation Assays indicate that test article, EDDN, did induce significant increases in the frequency of revertants for TA100 at 5000 µg/plate in the presence and absence of induced rat liver S-9 plus cofactors when compared to the solvent controls. This test followed OECD Test Guideline 471. Because the positive response was only observed at the highest dose, under the conditions of this test and according to the criteria set for evaluating the test results, the test article was weak positive in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay both with and without activation.

## TABLE OF CONTENTS

<b>STUDY DIRECTOR'S COMPLIANCE STATEMENT .....</b>	<b>2</b>
<b>QUALITY ASSURANCE UNIT'S STATEMENT .....</b>	<b>4</b>
<b>STUDY DIRECTOR'S SIGNATURE PAGE.....</b>	<b>5</b>
<b>ABSTRACT .....</b>	<b>6</b>
<b>INTRODUCTION .....</b>	<b>9</b>
<b>MATERIALS .....</b>	<b>11</b>
TEST ARTICLE	
CONTROL ARTICLES	
INDICATOR CELLS	
CULTURE CONDITIONS	
METABOLIC ACTIVATION SYSTEM	
<b>EXPERIMENTAL PROCEDURES .....</b>	<b>14</b>
DOCUMENTATION	
TEST SYSTEM IDENTIFICATION	
SOLUBILITY TEST	
PREPARATION OF TEST CULTURES	
PREPARATION OF METABOLIC ACTIVATION SYSTEM	
PREPARATION OF TEST ARTICLE DOSING SOLUTIONS	
RANGE FINDING TEST	
MUTATION ASSAYS	
<b>CRITERIA FOR A VALID ASSAY .....</b>	<b>21</b>
<b>EVALUATION OF TEST RESULTS .....</b>	<b>22</b>
<b>ARCHIVES.....</b>	<b>22</b>
<b>RESULTS .....</b>	<b>23</b>
SOLUBILITY TEST	
RANGE FINDING TEST	
MUTATION ASSAYS	
ANALYSIS OF DOSING SOLUTIONS	
<b>CONCLUSIONS .....</b>	<b>25</b>
<b>REFERENCES .....</b>	<b>25</b>
<b>APPENDIX I .....</b>	<b>26</b>
DATA TABLES	
<b>APPENDIX II.....</b>	<b>33</b>
DETAILED PLATE COUNTS AND BACKGROUND LAWN EVALUATION	



<b>APPENDIX III .....</b>	<b>46</b>
SITEK'S HISTORICAL POSITIVE AND SOLVENT CONTROL DATA	
<b>APPENDIX IV .....</b>	<b>50</b>
STUDY PROTOCOL, PROTOCOL AMENDMENTS AND PROTOCOL DEVIATION	
<b>APPENDIX V .....</b>	<b>71</b>
S-9 BATCH INFORMATION	

## INTRODUCTION

This study was conducted by Jian Song, Ph.D., Shashi Sharma, B.S. and Melkie Lulie, M.S. from August 11, 2009 to August 31, 2009, at SITEK Research Laboratories. The experimental procedures used to perform this study were essentially those of B. N. Ames, et al. (1), D. Maron and B. N. Ames (2), M. H. L. Green and W. J. Muriel (3), and S. Venitt and J. M. Parry (eds.) (4).

The purpose of this study was to evaluate the test article, Ethylenediamine dinitrate (EDDN), for its potential to cause mutations in the histidine operon of *Salmonella typhimurium* strains TA98, TA100, TA1535 and TA1537 and the tryptophan operon of *Escherichia coli* strain WP2 uvrA using the Plate Incorporation method of treatment. The *Salmonella typhimurium/Escherichia coli* Plate Incorporation Mutation Assay has been used extensively and has been demonstrated to be effective in detecting mutations caused by compounds from a wide range of chemical classes (1-4). Over several years, a large database of results has been accumulated which has confirmed its ability to detect genetically active compounds of most chemical classes with high efficiency (5).

The Ames Assay, in general, is performed using either the Plate Incorporation method or Pre-Incubation method. From the regulatory point of view, both assays are equally acceptable. In the Plate Incorporation method, treatment is performed by adding either 500  $\mu$ L of sterile, deionized water or 500  $\mu$ L of S-9 cofactor mix to tubes containing 2.0 mL of top agar supplemented with 1X histidine-biotin or 1X tryptophan solution. Immediately thereafter, 100  $\mu$ L of bacteria is added followed by 100  $\mu$ L of the appropriate test article concentration or solvent control. Each tube is vortexed for 2-3 seconds, and the contents are evenly distributed over a Vogel-Bonner bottom agar plate. Each plate is placed on a level surface until the top agar solidifies. The plates are inverted and incubated at  $37 \pm 1^\circ\text{C}$  for 48 to 72 hours. In the Pre-Incubation method the treatment is performed by adding either 500  $\mu$ L of sterile, deionized water or 500  $\mu$ L of S-9 cofactor mix to tubes followed by 100  $\mu$ L of bacteria and 100  $\mu$ L of the appropriate test article concentration or solvent. The tubes are incubated at  $37 \pm 1^\circ\text{C}$  for 20-30 minutes in a shaker incubator. Finally, 2.0 mL of top agar supplemented with 1X histidine-biotin or 1X tryptophan solution is added to the tube, the contents are vortexed 2-3 seconds and spread over a Vogel-Bonner bottom agar plate. The plates are inverted and incubated at  $37 \pm 1^\circ\text{C}$  for 48 to 72 hours. For some of the Azo and nitrosamine family of compounds, pre-incubation of the culture is required in order to be metabolized prior to plating. The Pre-Incubation method of treatment is performed at the request of the Sponsor.

The agar contains a trace of histidine that allows all the bacteria to undergo several divisions, thus producing a faint background lawn of bacteria. DNA replication is necessary in many cases for mutagenesis to occur and therefore the background lawn provides a good indicator of the inhibition of growth caused by the test chemical. Mutational events are rare, therefore it is essential that large populations of bacteria are used in mutagenicity testing. Maximum sensitivity is achieved by plating around  $1 \times 10^8$  bacteria.

The Ames Assay is the most widely used of all methods for determining the mutagenicity of chemicals. Because the bacterial strains used in this assay lack the enzymes necessary for metabolizing promutagens to ultimate mutagens, rat liver S-9 induced with Aroclor 1254 was added as a substitute for mammalian metabolism. This assay detects point mutations only and measures reverse mutation from acid auxotrophy to prototrophy. In this method, the bacterial strains used carry base substitution or frame shift mutations in operons coding for synthesis of specific amino acids. Therefore, these mutants (unlike their wild-type counterparts) cannot synthesize all their required amino acids from inorganic sources of nitrogen, being auxotrophic for the specific amino acids histidine and tryptophan. This assay determines whether the test article can reverse the effect of the pre-existing mutation by introducing a second mutation. When the cultures are exposed to a mutagen, some of the bacteria undergo genetic changes due to chemical interactions resulting in reversion of the bacteria to a non-histidine-requiring state or non-tryptophan-requiring state. The reverted bacteria will then grow in the absence of exogenous histidine or tryptophan thus providing an indication of the potential of the test chemical to cause mutation. Multiple tester strains are necessary because different strains are mutated by a different class (or different classes) of compound. The genotypes of the strains are verified concurrently.

The following are the details of possible mutations in the different strains (4):

Bacterial Strain	Mutation	Rfa	UvrB	R Factor (pKM101)	Type of Mutation
TA98*	HISD 3052	Yes	Yes	Yes	Frame shift
TA100**	HIS G46	Yes	Yes	Yes	Base Pair Substitution Frame shift
TA1535**	HIS G46 B-P	Yes	Yes	No	Base Pair Substitution
TA1537	HISC 3076	Yes	Yes	No	Frame shift
E. coli	Trp-	Yes	No (uvrA)	No	Base Pair Substitution

\* TA98 was derived from TA1538 (pKM101 plasmid added).

\*\* TA100 was derived from TA1535 (pKM101 plasmid added).

rfa - Defective lipopolysaccharide coat. More permeable to chemicals. (Sensitive to crystal violet.)

uvrB - Reduced error-free repair of some types of DNA damage. (Sensitive to UV light.)

R Factor (pKM101) - Increases sensitivity by enhancing error-prone DNA repair. (Ampicillin resistant if plasmid present.)

uvrA - Less DNA repair.

## MATERIALS

### TEST ARTICLE

1. Name:	<u>Ethylenediamine dinitrate (EDDN)</u>
2. CAS No.:	<u>Not Available</u>
3. Provided by:	<u>USA RDECOM, AMSRD-MSF</u> <u>Environmental Acquisition &amp; Logistics Sustaining</u> <u>Program</u> <u>Aberdeen Proving Ground, MD 21010</u>
4. Batch/Lot No.:	<u>ABY08L031S010</u>
5. Physical Description:	<u>Liquid(500mg/mL in deionized, distilled water)</u>
6. Shipping Conditions:	<u>Room Temperature</u>
7. Date Received at SITEK:	<u>August 5, 2009</u>
8. Storage Conditions:	<u>Refrigerated (1-5 °C)</u>
9. Purity:	<u>99.5%</u>
10. Expiration Date:	<u>Not Available</u>

The Certificate of Analysis is not available but the sponsor's data indicated the purity of the test article is 99.5%.

### CONTROL ARTICLES

#### Positive Controls

The positive control chemicals used for the tester strains in the presence and absence of exogenous metabolic activation are presented below:

<u>Strain</u>	<u>S-9</u>	<u>Chemical</u>	<u>Concentration (µg/plate)</u>
TA98	-	2-NF (2-Nitrofluorene)	5.0
TA98	+	2-AA (2-Aminoanthracene)	2.5
TA100	-	NaN <sub>3</sub> (Sodium Azide)	1.0
TA100	+	2-AA (2-Aminoanthracene)	2.5
TA1535	-	NaN <sub>3</sub> (Sodium Azide)	1.0
TA1535	+	2-AA (2-Aminoanthracene)	2.5
TA1537	-	9-AA (9-Aminoacridine)	75
TA1537	+	2-AA (2-Aminoanthracene)	5.0
WP2 uvrA	-	MMS (Methyl Methanesulfonate)	4000
WP2 uvrA	+	2-AA (2-Aminoanthracene)	20

The following is the information for each of the positive controls used in this assay:

<u>Chemical</u>	<u>*Source</u>	<u>CAS No.</u>	<u>Lot No.</u>	<u>Storage Conditions</u>	<u>Expiration Date</u>
2-AA	Aldrich	613-13-8	03403ED	1-5°C	08-04-14
9-AA	Aldrich	52417-22-8	1126KD	1-5°C	10-24-11
2-NF	Aldrich	607-57-8	092138A	1-5°C	03-23-12
NaAz	Sigma	26628-22-8	073K0119	1-5°C	03-23-12
MMS	Aldrich	66-27-3	06823KH	1-5°C	06-02-13

\* SIGMA-ALDRICH, St. Louis, MO 63178.

The positive controls 2-AA, 9-AA, and 2-NF were dissolved in DMSO. NaAz and MMS were dissolved in sterile deionized distilled water. Multiple vials of the above mentioned positive controls were prepared and frozen at  $-70^{\circ}\text{C} \pm 10$  were used in this assay. The source, lot number and expiration date of the DMSO used to prepare the positive controls are presented below:

Source: Sigma Chemical Company  
St. Louis, MO 63178

Lot No.: 10585CH

Storage Conditions: Room Temperature

Expiration Date: January 31, 2012

CAS No: 67-68-5

The source, batch numbers and expiration dates of the sterile deionized, distilled water (DD water) are presented below:

Source: SITEK Research Laboratories

Batch No.: 106, 107

Storage Conditions: Room Temperature

Expiration Dates: January 23, 2010 (Batch: 106), January 31, 2010 (Batch: 107)

### **Solvent Control**

The test article, EDDN, was prepared and diluted in DD water. Therefore, DD water was used as the solvent control. The source, batch numbers and expiration dates of the DD water are provided above.

## INDICATOR CELLS

### Source

The *Salmonella typhimurium* strains TA98, TA100, TA1535 and TA1537 were originally obtained from Dr. Bruce N. Ames, University of California, Berkeley. The *Escherichia coli* strain WP2 uvrA was obtained from Ms. Judy Mayo of Pharmacia and Upjohn Co., Kalamazoo, Michigan.

## CULTURE CONDITIONS

The cells were grown in Oxoid Nutrient Broth No. 2 (Oxoid LTD, Hampshire, England) in a shaker incubator rotating at approximately 120 rpm and maintained at a temperature of  $37 \pm 1^\circ\text{C}$ . Stock cultures of the tester strains were cryopreserved at SITEK Research Laboratories. Scrapes from the cryopreserved stock were used to initiate the overnight cultures for the test.

## METABOLIC ACTIVATION SYSTEM

For the activated portion of the range finding and mutation assays, the cells were exposed to the test article in conjunction with an exogenous metabolic activation system consisting of Aroclor-induced rat liver S-9 in 0.154M KCl plus cofactors (S-9 mix). The components of the standard S-9 mix were 8mM  $\text{MgCl}_2$ , 33mM KCl, 5mM glucose-6-phosphate, 4mM NADP, 100mM sodium phosphate buffer (pH 7.4), and 10% rat liver homogenate prepared from Aroclor 1254-induced, Sprague-Dawley rats. The S-9 batches used in this study were also evaluated for sterility, protein content and promutagen activity.

<u>Source:</u>	Molecular Toxicology, Inc., Boone, NC 28607
<u>Inducing Agent:</u>	Aroclor 1254
<u>S-9 Lot No.:</u>	2342
<u>Protein Content:</u>	35.5 mg/mL
<u>Storage Conditions:</u>	$\leq -70^\circ\text{C}$
<u>Expiration Date:</u>	September 9, 2010

<u>Source:</u>	AP Sciences, Inc., Columbia, MD 21045
<u>Inducing Agent:</u>	Phenobarbital/5,6-Benzoflavone
<u>S-9 Lot No.:</u>	01102-2001
<u>Protein Content:</u>	38.5 mg/mL
<u>Storage Conditions:</u>	$\leq -70^\circ\text{C}$
<u>Expiration Date:</u>	June 19, 2011

Detailed information concerning the S-9 batches used in the Assay is provided in Appendix V.

## **EXPERIMENTAL PROCEDURES**

### **DOCUMENTATION**

The materials, experimental procedures used in the performance of the study, experimental results and methods used in the evaluation of the results were documented in the study workbook.

### **TEST SYSTEM IDENTIFICATION**

#### **Plate Incorporation Method**

The Plate Incorporation method is performed by adding either 500 µL of sterile deionized, water or 500 µL of S-9 cofactor mix to tubes containing 2.0 mL of top agar supplemented with 1X histidine-biotin or 1X tryptophan solution. Immediately thereafter, 100 µL of respective bacteria is added followed by 100 µL of the appropriate test article concentration or solvent. Each tube is vortexed for 2-3 seconds, and the contents are evenly distributed over a Vogel-Bonner bottom agar plate. Each plate is placed on a level surface until the top agar solidified. The plates are inverted and incubated at  $37 \pm 1^{\circ}\text{C}$  for 48 to 72 hours.

#### **Labeling Plates for the Mutation Assay**

A sufficient number of Vogel-Bonner agar plates was removed from refrigerated storage and allowed to warm to room temperature. Each plate was then labeled with the following information: SITEK's test article number, experiment phase, presence or absence of rat liver S-9 mixture, concentration level code, and strain code. The following strain and concentration level codes were used:

##### **Strain Codes:**

1 = TA98      3 = TA1535      5 = WP2 uvrA  
2 = TA100      4 = TA1537

##### **Concentration Level Codes:**

0 = Solvent for the Test Article  
1 = 1st or highest Test Article concentration level  
2 = 2nd Test Article concentration level  
3 = 3rd Test Article concentration level  
4 = 4th Test Article concentration level  
5 = 5th Test Article concentration level or lowest Test Article concentration level for the  
Mutation Assays  
6 = 6th Test Article concentration level  
7 = 7th Test Article concentration level or lowest Test Article concentration level for the  
Range Finding Test.

In addition to the above, Mutation Assay viability plates that contained 10X histidine-biotin or 10X tryptophan were designated with the prefix "T".

#### **Labeling Positive Control Plates**

Vogel-Bonner agar plates were removed from refrigerated storage and allowed to warm to room temperature. Triplicate sets were labeled with the test article number, identity and concentration of the particular positive control, experimental phase, strain code, and the presence or absence of rat exogenous metabolic activation.

#### **Labeling Tester Strain Titer Plates**

Each tester strain titer plate was labeled with the following information: SITEK test article number, tester strain identity, and experimental phase and the prefix T.

#### **Labeling Tester Strain Characterization Plates**

##### **Histidine Requirement**

A single histidine-biotin plate was divided into four zones by drawing horizontal lines on the bottom of the plate with a marking pen and labeling each zone with a different *Salmonella* tester strain. A biotin-only control plate was labeled in a similar manner.

##### **rfa Mutation**

Nutrient agar plates were labeled with the *Salmonella* tester strain identification and "CV" (crystal violet).

##### **R-Factor**

A single ampicillin agar plate was labeled in a similar manner as the histidine-biotin plate.

##### **Tryptophan Requirement**

A tryptophan plate and a Vogel-Bonner agar control plate were labeled with the code for strain WP2 uvrA and used for confirmation of the tryptophan requirement.

#### **SOLUBILITY TEST**

The test article is liquid and the concentration is 500 mg/mL in deionized, distilled water. The sponsor specified water as the solvent of choice, therefore, a solubility test was not performed.



## PREPARATION OF TEST CULTURES

The methods used for the cryopreservation and cultivation of the tester strains are the procedures used by B. N. Ames et al. (1) as modified by D. Maron and B. N. Ames (2).

### Inoculation Procedures

Frozen ampules of strains TA98, TA100, TA1535, TA1537 and WP2 uvrA for the Mutation Assay were removed from liquid nitrogen and placed into crushed dry ice to prevent thawing. Scrapes were made using the tip of a sterile pipette, and these scrapes were transferred to a shaker flask containing approximately 50 mL of sterile Oxoid Nutrient Broth No. 2. The strains were incubated on a shaker at approximately 120 rpm and  $37 \pm 1^\circ\text{C}$ . The *Salmonella* strains were removed approximately 8-12 hours after the unit started and the *E. coli* strain was removed after approximately 4-6 hours.

### Harvesting Overnight Cultures

Before starting the experiment, the cultures were sampled and their percent transmittance (%T) was determined using a spectrophotometer set to a wavelength of 650 nm.

When the desired cell density of approximately  $5 \times 10^8$  to  $1 \times 10^9$  cells/mL (represented by a %T of between 25% and 10%, Optical Density of 0.6-1.0) was achieved, the cultures were placed on wet ice or kept at  $1-5^\circ\text{C}$  until needed.

## PREPARATION OF METABOLIC ACTIVATION SYSTEM

The S-9 cofactor mix was prepared as follows: For each mL of S-9 cofactor mix required, 0.335 mL of sterile deionized, distilled water was combined with 0.5 mL of 0.2M sodium phosphate buffer (pH 7.4), 0.04 mL of a 0.1M NADP solution, 5.0  $\mu\text{L}$  of 1M glucose-6-phosphate, and 0.02 mL of a 0.4M  $\text{MgCl}_2$ /1.65M KCl salt solution. This mixture was maintained on ice until just prior to use, whereupon 0.10 mL of S-9 in 0.154M KCl was added to the mixture.

## PREPARATION OF TEST ARTICLE DOSING SOLUTIONS

For the Range Finding Test and Definitive and Confirmatory Mutation Assays the test article was dissolved and diluted in the elected solvent in glass tubes. All the test article and control substance preparations and treatments were done under UV filtered lights to avoid possible problems of photoinactivation. The concentration and stability of the test article under experimental conditions was not determined.

## **RANGE FINDING TEST**

In order to determine the toxicity of the test article and to select appropriate test article concentrations for the Definitive Mutation Assay, a Range Finding Test was performed using strains TA100 and WP2 uvrA. The two strains have been successfully used and are sufficient to approximate the range of toxicity of the test article. Seven concentrations of the test article ranging from 5.0-5000 µg/plate were evaluated with and without induced rat liver S-9, using one plate per concentration.

### **Spontaneous Reversion Frequency**

Treatment was performed by adding either 500 µL of sterile deionized, distilled water or 500 µL of S-9 cofactor mix to tubes containing 2.0 mL of top agar supplemented with 1X histidine biotin or 1X tryptophan solution. Immediately thereafter, 100 µL of TA100 or WP2 uvrA was added followed by 100 µL of the appropriate test article concentration or solvent. Each tube was vortexed for 2-3 seconds, and the contents were evenly distributed over a Vogel-Bonner bottom agar plate. Each plate was placed on a level surface until the top agar solidified. The plates were inverted and incubated at  $37 \pm 1^{\circ}\text{C}$  for 48 to 72 hours.

### **Viable Count Determination**

Treatment and incubation were performed as described in the preceding paragraphs, except that approximately 250-500 cells of TA100 or WP2 uvrA were added to top agar supplemented with 10X histidine-biotin or 10X tryptophan solution.

After the incubation period was completed, the plates, starting with the highest test article concentration, were observed for the presence of precipitate. Plates having no interfering precipitate were counted for revertant colonies using an automatic colony counter (ARTEK Counter, Model 880, Manassas, Virginia 20110). Three counts were taken by rotating the plate on the counter stage and the median count was entered into a validated, MS Office Excel spreadsheet program designated as "2140A.xlw".

The background lawn was also evaluated. The following notations were used for the precipitate and background lawn evaluation:

#### **Chemical Precipitate:**

NP = No precipitate present.

SP = Slight precipitate - Noticeable compound on the plate; however, no influence on automated plate counting.

MP = Moderate precipitate - Moderate precipitate requiring hand counting for colony enumeration.

HP = Heavy precipitate – Large amount of compound on the plate rendering hand counting difficult.

#### **Background Lawn Evaluation:**

NL = Normal, healthy microcolony lawn.

SR = A noticeable thinning of the microcolony lawn compared to that of the solvent control plates.

MR = Marked thinning of the microcolony lawn and an increase in the size of the microcolonies compared to the solvent control plates.

ER = Extreme thinning of the microcolony lawn and a large increase in the size of the microcolonies compared to the solvent control plates.

AB = Absence of any microcolony bacterial lawn.

OP = Obscured by precipitate.

#### **Determination of Relative Cloning Efficiency**

The corrected viability counts from each concentration with and without activation in *Salmonella* strain TA100 and in *Escherichia coli* strain WP2 uvrA were compared with the respective solvent control viability counts. The resulting ratio is the Relative Cloning Efficiency (RCE) and was converted into a percentage, and the data were included in the Range Finding Test results. Relative Colony Efficiency measures the toxicity of test article in terms of cell viability. Generally, diluted cultures are treated at various test article concentrations and mixed with top agar containing higher concentration of respective amino acids (10X histidine-biotin or tryptophan). All viable bacteria are able to make countable colony. It is desirable, if possible, to test one or two higher concentrations around 50% toxicity level (reduction of RCE by 50% in comparison to concurrent control) in the mutation assays. This is not valid for the non-toxic test compound. Relative Colony Efficiencies are not determined during the Definitive and Confirmatory Mutation Assays as the range of toxicity information is already available from the Range Finding Assay.

### **MUTATION ASSAYS**

#### **Definitive Mutation Assay**

Concentrations for the Definitive Mutation Assay were selected based on the results of the Range Finding Test. The Definitive Mutation Assay was performed with the four *Salmonella typhimurium* tester strains (TA98, TA100, TA1535 and TA1537) and *Escherichia coli* strain WP2 uvrA using the plate incorporation method of treatment. Treatment was performed by adding either 500 µL of sterile deionized, distilled water or 500 µL of rat S-9 cofactor mix to tubes containing 2.0 mL of top agar supplemented with 1X histidine-biotin or 1X tryptophan solution. Immediately

thereafter, 100  $\mu$ L of strains TA98, TA100, TA1535, TA1537 or WP2 uvrA were added, followed by 100  $\mu$ L of the appropriate test article concentration or solvent. The positive controls were treated with 100  $\mu$ L of the appropriate stock solutions. Each tube was vortexed for 2-3 seconds and the contents were evenly distributed over a Vogel-Bonner bottom agar plate. Each plate was placed on a level surface until the top agar solidified. The plates then were inverted and incubated at  $37 \pm 1^\circ\text{C}$  for approximately 48 - 72 hours.

### **Tester Strain Titer Determination**

Each tester strain was diluted to determine the approximate number of viable cells delivered to the assay plates. Therefore, approximately 250-500 cells were added to top agar supplemented with 10X histidine-biotin or 10X tryptophan solution. Each tube was vortexed for 2-3 seconds and the contents were evenly distributed on bottom agar plates. The plates were incubated at  $37 \pm 1^\circ\text{C}$  for approximately 48 to 72 hours.

### **Tester Strain Characterization**

All of the *Salmonella typhimurium* strains used in the assay were confirmed for the histidine requirement and the rfa mutation. In addition, strains TA98 and TA100 were tested for the presence of the pKM101 plasmid. *Escherichia coli* strain WP2 uvrA was confirmed for the tryptophan requirement.

### **Histidine or Tryptophan Requirement**

A streak of each tester strain was made by dipping a sterilized micro pipette tip into the appropriate undiluted tester strain suspension and drawing it across the surface in the appropriate region of a labeled histidine-biotin or tryptophan plate, as well as control plates. The plates were incubated at  $37 \pm 1^\circ\text{C}$  for approximately 12 to 24 hours.

### **rfa Mutation**

For each of the *Salmonella* tester strains, a 100  $\mu$ L aliquot of the undiluted culture was added to a tube containing 2.0 mL of 1X histidine-biotin solution in top agar. Each tube was vortexed for 2-3 seconds, and the contents were poured onto an appropriately labeled nutrient agar plate. After allowing the plate to solidify, a sterile disc was aseptically placed in the center of the agar overlay. Ten  $\mu$ L of a 1.0 mg/mL crystal violet solution was then added to the disc. The plates were incubated at  $37 \pm 1^\circ\text{C}$  for approximately 12 to 24 hours.

### **R-Factor Plasmid**

For each of the *Salmonella* tester strains, a 100  $\mu$ L aliquot of the undiluted culture was added to a tube containing 2.0 mL of 1X histidine-biotin solution in top agar. Each tube was vortexed for 2-3 seconds, and the contents were poured onto an appropriately labeled nutrient agar plate. After allowing the plate to solidify, a sterile disc was aseptically placed in the center

of the agar overlay. Ten  $\mu\text{L}$  of a 0.8% ampicillin solution was then added to the disc. The plates were incubated at  $37 \pm 1^\circ\text{C}$  for approximately 12 to 24 hours.

### **uvrB Deletion**

After the cryopreservation of the *Salmonella typhimurium* strains and the *E. coli* strain, the stock ampules were checked for *uvrB* deletion. For each of the *Salmonella* tester strains, a 100  $\mu\text{L}$  aliquot of the undiluted culture was added to a tube containing 2.0 mL of 1X histidine-biotin solution top agar. One hundred  $\mu\text{L}$  of the *E. coli* strain was added to a tube containing 2.0 mL of 1X tryptophan solution top agar. Each tube was vortexed for 2-3 seconds, and the contents were poured onto an appropriately labeled nutrient agar plate. After allowing the plate to solidify, half of the plate was covered with foil. The plates were placed under UV light for thirty seconds and then incubated at  $37 \pm 1^\circ\text{C}$  for approximately 12 to 24 hours.

### **Evaluation of Assay Results**

After the incubation period was completed, the plates, starting with the highest test article concentration, were observed for the presence of precipitate. Plates were counted for the frequency of revertant colonies using an ARTEK counter, model 880. Three counts were taken by rotating the plate on the counter stage and the median count was entered into a validated, MS Office Excel spreadsheet program designated as "2140B.xlw".

The background lawn was also evaluated. The same notations as in the Range Finding Test were used to evaluate the precipitate and background lawn.

### **Evaluation of Tester Strain Characterization**

The requirement for histidine or tryptophan was demonstrated by the growth of the tester strains on plates supplemented with histidine or tryptophan and the lack of growth on the control plates.

The presence of the *rfa* mutation was evaluated by measuring the zone of inhibition around the crystal violet disc. A zone  $\geq 12$  mm in diameter was evidence of appropriate inhibition.

The presence of the pKM101 plasmid was demonstrated by the growth of strains TA98 and TA100 and the lack of growth of strains TA1535 and TA1537 by observation the zone of inhibition around the ampicillin disc.

### **Tabulation of Colony Counts**

The colony counts provided by the automatic colony counter or by hand count were raw counts and were not corrected to reflect actual counts. Correction of the counts was performed by computer. The data tables presented in Appendix I contain the corrected values. The correction factor was determined by comparing a wide range of manual and automatic counts, as described

in SITEK's SOP No. 21.0. The relationship was linear, and the counts were corrected by using the following formula:

$$\text{Corrected Count} = (\text{Raw Counts}) (1.0571607) + 3.09496$$

### **Confirmatory Mutation Assay**

If the first Mutation Assay gives negative or equivocal results, a confirmatory Mutation Assay will be performed. The test article treatment concentrations may be altered based on the results obtained in the first Mutation Assay. On the other hand, if the results of the first Mutation Assay are clearly positive, a confirmatory Mutation Assay may or may not be performed depending on the Sponsor's instructions.

## **CRITERIA FOR A VALID ASSAY**

The following criteria were used as guidelines in evaluating the acceptability of the Mutation Assay. Because it is impossible to formulate criteria that would apply to every configuration of data generated by the assay, the Study Director was responsible for the ultimate decision regarding the acceptability of the results.

### **Solvent Control Cultures**

The mean reversion frequency (number of colonies on Agar plates) of the test article solvent control plates for each tester strain should fall within the following ranges:

TA98	30 ± 15	WP2 uvrA	15 ± 10
TA100	100 ± 70		
TA1535	20 ± 15		
TA1537	15 ± 12		

### **Positive Controls**

The results for the positive control cultures were considered acceptable if the treated strains had a mean reversion frequency that was three times or higher, than the mean reversion frequency of the solvent control plates.

### **Tester Strain Characterization**

All of the *Salmonella typhimurium* strains were confirmed positive for histidine dependence. *Escherichia coli* strain WP2 uvrA was confirmed positive for tryptophan dependence.

All of the *Salmonella typhimurium* strains were confirmed positive for the rfa mutation as evidenced by sensitivity to crystal violet.

The R-factor strains, TA98 and TA100, were confirmed positive for the pKM101 plasmid

as evidenced by ampicillin resistance.

The titer of the stock cultures for each strain indicated that the stock cultures contained greater than  $0.5 \times 10^9$  bacteria per mL.

## EVALUATION OF TEST RESULTS

The following criteria were used as guidelines in evaluating the results of the Mutation Assay for a negative, positive or equivocal response. Because it is impossible to write criteria that would apply to every configuration of data generated by the assay, the Study Director was responsible for the ultimate decision concerning the results.

### **Criteria for a Negative Response**

A response was considered to be negative if all of the strains treated with the test article had mean reversion frequencies that were less than twice that of the mean reversion frequencies of the corresponding solvent control plates in TA98 and TA100 and less than three times in TA1535, TA1537 and WP2 uvrA, and there was no evidence of a concentration-dependent response.

### **Criteria for a Positive Response**

A response was considered to be positive if either strain TA98 or TA100 exhibited a mean reversion frequency that was at least double the mean reversion frequency of the corresponding solvent control in at least one concentration, or if either strain TA1535, TA1537 or WP2 uvrA exhibited a three-fold increase in the mean reversion frequency compared to the solvent control in at least one concentration. In addition, the response must have been concentration-dependent or increasing concentrations of the test article must have showed increasing mean reversion frequencies. In evaluating the results, consideration was given to the degree of toxicity exhibited by the concentration causing the 2 to 3-fold or greater increase in reversion frequency and the magnitude of the increase in reversion frequency.

### **Criteria for an Equivocal Response**

A response was considered equivocal if it did not fulfill the criteria of either a negative or a positive response and/or the Study Director did not consider the response to be either positive or negative.

## ARCHIVES

The raw data, documentation, protocol, protocol amendment/deviation and a copy of the Final Report, along with an electronic file containing data tables and the Final Report of the study, will be maintained by SITEK Research Laboratories at Dr. Kirby's private residence until arrangements are made to transfer them to the Sponsor.

## RESULTS

### SOLUBILITY TEST

Solubility test was not performed because Ethylenediamine dinitrate (EDDN) for this assay is 500 mg/mL in deionized, distilled water and the solvent is deionized, distilled water.

### RANGE FINDING TEST

Summaries of the results of the Range Finding Test are presented in Tables 1 and 2 (Appendix I). The individual plate counts and background lawn evaluations are presented in Appendix II.

#### TA100:

The Relative Cloning Efficiencies (RCEs) at the concentrations of 5.0 to 5000 µg/plate without activation ranged from 63% to 134%. In the presence of the activation system, the RCEs at the concentrations of 5.0 to 5000 µg/plate ranged from 145% to 270%. No significant decrease in RCE or the number of revertants or thinning of background lawn was observed at all dose levels. The test article was nontoxic at all concentrations tested. No precipitate was observed at any of the test concentrations.

#### WP2 uvrA:

The Relative Cloning Efficiencies (RCEs) at the concentrations of 5.0 to 5000 µg/plate without activation ranged from 72% to 117%. In the presence of the activation system, the RCEs at the concentrations of 5.0 to 5000 µg/plate ranged from 76% to 119%. No significant decrease in RCE or the number of revertants was observed at all dose levels. The test article was nontoxic at all concentrations tested. No precipitate or thinning of background lawn was observed at any of the test concentrations.

### MUTATION ASSAYS

#### Definitive Mutation Assay

Summaries of the results of the Definitive Mutation Assay are presented in Tables 3 and 4 in Appendix I. The individual plate counts and background lawn evaluations are presented in Appendix II.

The Definitive Mutation Assay (B1), using the plate incorporation method, 100, 500, 1000, 3000 and 5000 µg/plate were test with the four *Salmonella* tester strains (TA98, TA100, TA1535 and TA1537) and with *E. coli* strain WP2 uvrA. The results indicate that under the condition of this study, EDDN, the results indicate that, EDDN, did cause positive increases in the mean number of revertants per plates with the tester strain TA100 without activation (3.1-fold) at 5000 µg/plate. But no significant dose-dependent response was observed. No positive increases were observed with other tester strain/activation combinations. The background lawns were normal. The solvent controls and positive controls fulfilled the requirements of a valid test. The results



indicate the test article was equivocal response in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay without activation.

### **Confirmatory Mutation Assay**

Summaries of the results of the Confirmatory Mutation Assay are presented in Tables 5 and 6 in Appendix I. The individual plate counts and background lawn evaluations are presented in Appendix II.

The Confirmatory Mutation Assay (B2) was performed using the plate incorporation method at concentrations of 3000, 3500, 4000, 4500 and 5000 µg/plate for *Salmonella typhimurium* (TA98, TA100, TA1535 and TA1537) and *Escherichia coli* both with and without activation because a positive response was only observed at 5000 µg/plate for TA100 without activation in the Definitive Assay. The results indicate that, EDDN, did cause positive increased in the mean number of revertants per plates with the tester strain TA100 without activation (2.6-fold) at 5000 µg/plate, also with activation was positive (2.5-fold) at 5000 µg/plate in the Confirmatory Mutation Assay and dose-dependent response was observed. No positive increases were observed with other tester strain/activation combinations. The background lawns were normal. The solvent controls and positive controls fulfilled the requirements of a valid test. The results indicate the test article was weak positive in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay both with and without activation..

SITEK's historical data for positive and solvent controls are presented in Appendix III.

### **ANALYSIS OF DOSING SOLUTIONS**

The Sponsor did not elect to have dosing solutions analyzed.

## CONCLUSIONS

The test article, Ethylenediamine dinitrate (EDDN, 99.5% pure) was tested in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay in the presence and absence of induced rat liver S-9. Definitive and Confirmatory Assays were performed.

The results of the Mutation Assays indicate that test article, EDDN, induced significant increases in the revertant frequencies for the tester strains TA100 without activation in the Definitive Assay and both without and with activation in the Confirmatory Assay only at the highest concentration which is 5000 µg/plate when compared to the solvent controls. Dose-dependent response was only observed in the Confirmatory Assay.

Therefore, under the conditions of this study, the test article, EDDN, was weak positive in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay both with and without activation. The strength and stability of the test article, dosing solutions, under the experimental conditions, were not determined and the impact on results and the conclusion is not known.

## REFERENCES

1. Ames, B. N., J. McCann and E. Yamasaki. Methods for detecting carcinogens and mutagens with the *Salmonella*/mammalian-microsome mutagenicity test. *Mut. Res.*, 31:347-367, 1975.
2. Maron, D., and B. N. Ames. Revised methods for the *Salmonella* mutagenicity test. *Mut. Res.*, 113:173-215, 1983.
3. Green, M. H. L., and W. J. Muriel. Mutagen testing using trp<sup>+</sup> reversion in *Escherichia coli*. In: B. J. Kilbey, et al. (eds.), *Handbook of Mutagenicity Test Procedures*, pp. 65-94, Elsevier North Holland Biomedical Press, Amsterdam, 1977.
4. Venitt, S., and J. M. Parry (eds.). *Mutagenicity Testing: A Practical Approach*. IRL Press, Oxford, England and Washington, D.C., 1984.
5. Jagannath, D.R. & Brusick, D.J. (1983) Mutagens and Carcinogens in Bacteria. In *A Guide to General Toxicology* (eds. Homburger, F., Hayes, J.A. & Pelikan, E.W.), S. Karger AG, P.O. Box, CH-4009 Basel, Switzerland, pp 344-364.

**APPENDIX I**  
**DATA TABLES**

TABLE 1

**SALMONELLA TYPHIMURIUM PLATE INCORPORATION MUTATION ASSAY**  
**RANGE FINDING TEST RESULTS**

SPONSOR: USA RDECOM, AMSRD-MSF      SITEK STUDY NO.: 1003-2140  
 EXPERIMENT NO.: A1      SOLVENT: DD Water  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)      STRAIN: TA100

WITHOUT ACTIVATION						WITH S-9 ACTIVATION					
Test Article Conc. µg/Plate	No. of Rever- tants/ Plate	Chem. PPT. Eval.*	Back-ground Lawn Evaluation**	No. of Viable Colo- nies/ Plate	Rela- tive Cloning Effi- ciency (RCE)	Test Article Conc. µg/Plate	No. of Rever- tants/ Plate	Chem. PPT. Eval.*	Back-ground Lawn Evaluation**	No. of Viable Colo- nies/ Plate	Rela- tive Cloning Effi- ciency (RCE)
5.0	111	NP	NL	395	124%	5.0	61	NP	NL	225	145%
10	60	NP	NL	428	134%	10	51	NP	NL	238	154%
50	100	NP	NL	390	122%	50	76	NP	NL	243	157%
100	101	NP	NL	406	127%	100	61	NP	NL	235	152%
500	87	NP	NL	299	94%	500	70	NP	NL	268	173%
1000	108	NP	NL	280	88%	1000	88	NP	NL	341	220%
5000	172	NP	NL	201	63%	5000	108	NP	NL	419	270%
SOLV. CONT.	100	NP	NL	319	100%	SOLV. CONT.	76	NP	NL	155	100%

$$\text{RCE} = \frac{\text{No. of Colonies in Test Plates}}{\text{No. of Colonies in Solvent Control Plates}} \times 100$$

## \* Chemical Precipitate Evaluation

NP = No precipitate

SP = Slight precipitate; noticeable precipitate on the plate, but no interference with automated plate counting

MP = Moderate precipitate; marked precipitate necessitating hand counting for colony enumeration

HP = Heavy precipitate; large amount of precipitate rendering hand counting difficult or impossible

## \*\* Background Lawn Evaluation

NL = Normal, healthy microcolony lawn

SR = Noticeable thinning of the microcolony lawn compared to control

MR = Marked thinning of the microcolony lawn and increase in size of microcolonies compared to control

ER = Extreme thinning of the microcolony lawn and large increase in size of microcolonies compared to control

AB = Absence of microcolonies

OP = Obscured by precipitate

Verified by: QA PK SD 2/25/10

TABLE 2

**ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**RANGE FINDING TEST RESULTS**

SPONSOR: USA RDECOM, AMSRD-MSF      SITEK STUDY NO.: 1003-2140  
 EXPERIMENT NO.: A1      SOLVENT: DD Water  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)      STRAIN: WP2 uvrA

WITHOUT ACTIVATION						WITH S-9 ACTIVATION					
Test Article Conc. µg/Plate	No. of Rever- tants/ Plate	Chem. PPT. Eval.*	Back-ground Lawn Evalu- ation**	No. of Viable Colo- nies/ Plate	Rela- tive Cloning Effi- ciency (RCE)	Test Article Conc. µg/Plate	No. of Rever- tants/ Plate	Chem. PPT. Eval.*	Back-ground Lawn Evalu- ation**	No. of Viable Colo- nies/ Plate	Rela- tive Cloning Effi- ciency (RCE)
5.0	13	NP	NL	1074	117%	5.0	14	NP	NL	1293	119%
10	13	NP	NL	737	80%	10	8	NP	NL	826	76%
50	5	NP	NL	759	82%	50	20	NP	NL	1063	97%
100	9	NP	NL	729	79%	100	9	NP	NL	881	81%
500	15	NP	NL	678	74%	500	16	NP	NL	905	83%
1000	9	NP	NL	993	108%	1000	13	NP	NL	1103	101%
5000	10	NP	NL	663	72%	5000	13	NP	NL	867	79%
SOLV. CONT.	12	NP	NL	921	100%	SOLV. CONT.	20	NP	NL	1091	100%

$$\text{RCE} = \frac{\text{No. of Colonies in Test Plates}}{\text{No. of Colonies in Solvent Control Plates}} \times 100$$

## \* Chemical Precipitate Evaluation

NP = No precipitate

SP = Slight precipitate; noticeable precipitate on the plate, but no interference with automated plate counting

MP = Moderate precipitate; marked precipitate necessitating hand counting for colony enumeration

HP = Heavy precipitate; large amount of precipitate rendering hand counting difficult or impossible

## \*\* Background Lawn Evaluation

NL = Normal, healthy microcolony lawn

SR = Noticeable thinning of the microcolony lawn compared to control

MR = Marked thinning of the microcolony lawn and increase in size of microcolonies compared to control

ER = Extreme thinning of the microcolony lawn and large increase in size of microcolonies compared to control

AB = Absence of microcolonies

OP = Obscured by precipitate

Verified by: QA PK SD 2/25/10

TABLE 3  
SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
 DEFINITIVE MUTATION ASSAY RESULTS - WITHOUT ACTIVATION

SPONSOR: USA RDECOM, AMSRD-MSF  
 EXPERIMENT NO.: B-1  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SITEK STUDY NO.: 1003-2140  
 SOLVENT: DD Water  
 CONC. IN: µg/plate

<u><i>S. typhimurium</i></u>		Average No. of Revertants Per Plate						
		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDDED: 9.440E+07	REVERTANTS	280	21	18	23	30	31	29
	STD. DEV.	52	2	1	6	3	4	7
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDDED: 8.380E+07	REVERTANTS	239	56	65	78	57	73	174
	STD. DEV.	21	5	11	3	10	8	87
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDDED: 1.824E+08	REVERTANTS	569	22	23	18	16	20	20
	STD. DEV.	31	5	1	3	8	2	3
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDDED: 1.418E+08	REVERTANTS	185	9	13	13	13	14	16
	STD. DEV.	136	1	3	3	3	2	4
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u><i>E. coli</i></u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDDED: 2.864E+08	REVERTANTS	379	17	16	14	14	16	13
	STD. DEV.	63	5	2	4	1	4	1
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

Verified by: QA RK SD 2/25/10

TABLE 4  
SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
 DEFINITIVE MUTATION ASSAY RESULTS - WITH S-9 ACTIVATION

SPONSOR: USA RDECOM, AMSRD-MSF  
 EXPERIMENT NO.: B-1  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SITEK STUDY NO.: 1003-2140  
 SOLVENT: DD Water  
 CONC. IN: µg/plate

<u>S. typhimurium</u>		Average No. of Revertants Per Plate						
		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDED: 9.440E+07	REVERTANTS	888	25	19	20	19	19	23
	STD. DEV.	115	3	3	3	3	4	4
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDED: 8.380E+07	REVERTANTS	967	78	59	87	81	109	128
	STD. DEV.	127	3	12	11	6	12	56
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.824E+08	REVERTANTS	217	16	19	16	15	12	20
	STD. DEV.	11	1	2	3	2	6	6
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.418E+08	REVERTANTS	166	12	10	9	14	13	15
	STD. DEV.	8	3	3	2	8	1	1
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	REVERTANTS	498	17	18	16	18	23	26
	STD. DEV.	26	4	2	4	4	9	5
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

Verified by: QA RK SD 2/25/10

TABLE 5  
SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
 CONFIRMATORY MUTATION ASSAY RESULTS - WITHOUT ACTIVATION

SPONSOR: USA RDECOM, AMSRD-MSF  
 EXPERIMENT NO.: B-2  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SITEK STUDY NO.: 1003-2140  
 SOLVENT: DD Water  
 CONC. IN: µg/plate

<i>S. typhimurium</i>		Average No. of Revertants Per Plate						
		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: TA98 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.078E+08	REVERTANTS	357	20	26	29	34	50	51
	STD. DEV.	93	6	8	5	6	14	6
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA100 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.460E+08	REVERTANTS	251	82	51	65	92	121	213
	STD. DEV.	21	16	4	9	11	19	63
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1535 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.420E+08	REVERTANTS	535	21	19	19	18	25	25
	STD. DEV.	32	2	3	4	5	9	7
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDED: 3.088E+08	REVERTANTS	213	12	7	7	8	8	9
	STD. DEV.	26	5	1	2	1	1	3
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<i>E. coli</i>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: WP2 <i>uvrA</i> DATE PLATED: 08/27/2009 CELLS SEEDED: 2.022E+08	REVERTANTS	469	15	20	22	17	21	17
	STD. DEV.	34	3	7	2	4	0	5
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

Verified by: QA PK SD 2/25/10



TABLE 6  
SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
 CONFIRMATORY MUTATION ASSAY RESULTS - WITH S-9 ACTIVATION

SPONSOR: USA RDECOM, AMSRD-MSF      SITEK STUDY NO.: 1003-2140  
 EXPERIMENT NO.: B-2      SOLVENT: DD Water  
 TEST ARTICLE: Ethylenediamine dinitrate (EDDN)      CONC. IN: µg/plate

<u>S. typhimurium</u>		Average No. of Revertants Per Plate						
		Positive Control	Solvent Control	Concentration per plate				
		3000	3500	4000	4500	5000		
STRAIN: TA98 DATE PLATED: 08/27/2009  CELLS SEEDED: 2.078E+08	REVERTANTS	975	37	31	27	34	37	53
	STD. DEV.	70	8	6	2	3	1	8
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA100 DATE PLATED: 08/27/2009  CELLS SEEDED: 2.460E+08	REVERTANTS	965	97	110	74	100	115	243
	STD. DEV.	125	15	27	6	12	12	93
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1535 DATE PLATED: 08/27/2009  CELLS SEEDED: 2.420E+08	REVERTANTS	198	16	28	23	27	16	17
	STD. DEV.	39	1	5	7	6	5	3
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1537 DATE PLATED: 08/27/2009  CELLS SEEDED: 3.088E+08	REVERTANTS	250	11	10	12	11	12	8
	STD. DEV.	23	2	4	5	6	3	2
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
		3000	3500	4000	4500	5000		
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009  CELLS SEEDED: 2.022E+08	REVERTANTS	397	20	17	20	16	16	16
	STD. DEV.	17	5	1	2	1	4	3
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

Verified by: QA RK SD 2/25/10

**APPENDIX II**  
**DETAILED PLATE COUNTS AND**  
**BACKGROUND LAWN EVALUATION**

## SALMONELLA TYPHIMURIUM PLATE INCORPORATION MUTATION ASSAY

EXPERIMENT NO.:

A1

SITEK STUDY NO.: 1003-2140

**TEST ARTICLE:**

### Ethylenediamine dinitrate (EDDN)

**SOLVENT:** DD Water

STRAIN: TA100

WITHOUT ACTIVATION

Test Article Conc. $\mu\text{g}/\text{Plate}$	No. of Revertants Per Plate (raw) (corrected)		Chem. Background PPT. Lawn Eval.* Evaluation**		No. of Viable Colonies/Plate (raw) (corrected)		Relative Cloning Efficiency (RCE)
5.0	102	111	NP	NL	371	395	124%
10	54	60	NP	NL	402	428	134%
50	92	100	NP	NL	366	390	122%
100	93	101	NP	NL	381	406	127%
500	79	87	NP	NL	280	299	94%
1000	99	108	NP	NL	262	280	88%
5000	160	172	NP	NL	187	201	63%
SOLVENT CONTROL	92	100	NP	NL	299	319	100%

$$\text{RCE} = \frac{\text{No. of Colonies in Test Plates}}{\text{No. of Colonies in Solvent Control Plates}} \times 100$$

### \* Chemical Precipitate Evaluation

NP = No precipitate

SP = Slight precipitate; noticeable precipitate on the plate, but no interference with automated plate counting

MP = Moderate precipitate; marked precipitate necessitating hand counting for colony enumeration

HP = Heavy precipitate; large amount of precipitate rendering hand counting difficult or impossible

## \*\* Background Lawn Evaluation

NL = Normal, healthy microcolony lawn

SR = Noticeable thinning of the microcolony lawn compared to control

MR = Marked thinning of the microcolony lawn and increase in size of microcolonies compared to control

ER = Extreme thinning of the microcolony lawn and large increase in size of microcolonies compared to control

**AB = Absence of microcolonies**

OP = Obscured by precipitate

**SALMONELLA TYPHIMURIUM PLATE INCORPORATION MUTATION ASSAY**  
**RANGE FINDING TEST COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.:	A1	SITEK STUDY NO.:	1003-2140
TEST ARTICLE:	Ethylenediamine dinitrate (EDDN)	SOLVENT:	DD Water
		STRAIN:	TA100

**WITH S-9 ACTIVATION**

Test Article Conc. $\mu\text{g}/\text{Plate}$	No. of Revertants Per Plate (raw) (corrected)		Chem. Background PPT. Lawn Eval.* Evaluation**		No. of Viable Colonies/Plate (raw) (corrected)		Relative Cloning Efficiency (RCE)
5.0	55	61	NP	NL	210	225	145%
10	45	51	NP	NL	222	238	154%
50	69	76	NP	NL	227	243	157%
100	55	61	NP	NL	219	235	152%
500	63	70	NP	NL	251	268	173%
1000	80	88	NP	NL	320	341	220%
5000	99	108	NP	NL	393	419	270%
SOLVENT CONTROL	69	76	NP	NL	144	155	100%

$$\text{RCE} = \frac{\text{No. of Colonies in Test Plates}}{\text{No. of Colonies in Solvent Control Plates}} \times 100$$

### \* Chemical Precipitate Evaluation

NP = No precipitate

SP = Slight precipitate; noticeable precipitate on the plate, but no interference with automated plate counting

MP = Moderate precipitate; marked precipitate necessitating hand counting for colony enumeration

HP = Heavy precipitate; large amount of precipitate rendering hand counting difficult or impossible

## \*\* Background Lawn Evaluation

NL = Normal, healthy microcolony lawn

SR = Noticeable thinning of the microcolony lawn compared to control

MR = Marked thinning of the microcolony lawn and increase in size of microcolonies compared to control

ER = Extreme thinning of the microcolony lawn and large increase in size of microcolonies compared to control

**AB = Absence of microcolonies**

OP = Obscured by precipitate

### ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY

## RANGE FINDING TEST COLONY COUNTS AND BACKGROUND LAWN EVALUATION

EXPERIMENT NO.:	A1	SITEK STUDY NO.:	1003-2140
TEST ARTICLE:	Ethylenediamine dinitrate (EDDN)	SOLVENT:	DD Water
		STRAIN:	WP2 uvrA

### WITHOUT ACTIVATION

Test Article Conc. µg/Plate	No. of Revertants		Chem. Background		No. of Viable Colonies/Plate		Relative Cloning Efficiency (RCE)
	Per Plate		PPT.	Lawn			
	(raw)	(corrected)	Eval.*	Evaluation**	(raw)	(corrected)	
5.0	9	13	NP	NL	1013	1074	117%
10	9	13	NP	NL	694	737	80%
50	2	5	NP	NL	715	759	82%
100	6	9	NP	NL	687	729	79%
500	11	15	NP	NL	638	678	74%
1000	6	9	NP	NL	936	993	108%
5000	7	10	NP	NL	624	663	72%
SOLVENT CONTROL	8	12	NP	NL	868	921	100%

$$\text{RCE} = \frac{\text{No. of Colonies in Test Plates}}{\text{No. of Colonies in Solvent Control Plates}} \times 100$$

### \* Chemical Precipitate Evaluation

NP = No precipitate

SP = Slight precipitate; noticeable precipitate on the plate, but no interference with automated plate counting

MP = Moderate precipitate; marked precipitate necessitating hand counting for colony enumeration

HP = Heavy precipitate; large amount of precipitate rendering hand counting difficult or impossible

## \*\* Background Lawn Evaluation

NL = Normal, healthy microcolony lawn

SR = Noticeable thinning of the microcolony lawn compared to control

MR = Marked thinning of the microcolony lawn and increase in size of microcolonies compared to control

ER = Extreme thinning of the microcolony lawn and large increase in size of microcolonies compared to control

**AB = Absence of microcolonies**

OP = Obscured by precipitate



SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
DEFINITIVE MUTATION ASSAY RAW COLONY COUNTS AND BACKGROUND LAWN EVALUATION

EXPERIMENT NO.: B-1

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITHOUT ACTIVATION

<u>S. typhimurium</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDED: 9.440E+07	REVERTANTS	312	18	13	13	25	23	24
	PER	261	15	14	22	27	27	31
	PLATE	214	19	14	23	23	29	19
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDED: 8.380E+07	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	246	51	53	74	44	62	248
	PER	209	54	70	69	47	61	154
	PLATE	215	45	52	69	62	75	84
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.824E+08	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	518	14	18	12	19	17	13
	PER	519	17	20	18	5	17	19
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.418E+08	PLATE	569	23	20	13	12	14	17
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	92	6	11	12	11	11	15
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.418E+08	PER	104	5	11	7	10	8	8
	PLATE	321	6	7	9	6	12	13
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	REVERTANTS	372	12	14	6	11	16	8
	PER	406	8	12	13	10	11	8
	PLATE	289	18	11	11	9	9	10
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**DEFINITIVE MUTATION ASSAY RAW COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-1

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITH S-9 ACTIVATION

<u><i>S. typhimurium</i></u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDED: 9.440E+07	REVERTANTS	917	19	14	12	14	13	22
	PER	714	24	18	17	18	19	15
	PLATE	881	19	13	18	13	12	20
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDED: 8.380E+07	REVERTANTS	812	69	44	89	76	111	92
	PER	879	71	49	82	78	88	84
	PLATE	1045	74	65	68	68	101	180
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.824E+08	REVERTANTS	215	13	16	11	9	14	22
	PER	199	12	17	15	12	8	13
	PLATE	195	11	13	10	12	4	12
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.418E+08	REVERTANTS	156	11	4	8	6	9	11
	PER	146	7	9	5	19	10	11
	PLATE	160	7	7	5	6	9	10
	LAWN	NL	NL	NL	NL	NL	NL	NL
		PRECIPITATE	NP	NP	NP	NP	NP	NP

<u><i>E. coli</i></u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	REVERTANTS	459	8	16	10	12	12	17
	PER	496	14	14	9	18	28	26
	PLATE	449	16	12	16	12	16	21
	LAWN	NL	NL	NL	NL	NL	NL	NL
		PRECIPITATE	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.



**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**DEFINITIVE MUTATION ASSAY CORRECTED COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-1

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

## WITHOUT ACTIVATION

<u>S. typhimurium</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDED: 9.440E+07	REVERTANTS	333	22	17	17	30	27	28
	PER	279	19	18	26	32	32	36
	PLATE	229	23	18	27	27	34	23
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDED: 8.380E+07	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	263	57	59	81	50	69	265
	PER	224	60	77	76	53	68	166
	PLATE	230	51	58	76	69	82	92
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.824E+08	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	551	18	22	16	23	21	17
	PER	552	21	24	22	8	21	23
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDED: 1.418E+08	PLATE	605	27	24	17	16	18	21
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	100	9	15	16	15	15	19
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	PER	113	8	15	10	14	12	12
	PLATE	342	9	10	13	9	16	17
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	REVERTANTS	396	16	18	9	15	20	12
	PER	432	12	16	17	14	15	12
	PLATE	309	22	15	15	13	13	14
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDED: 2.864E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**DEFINITIVE MUTATION ASSAY CORRECTED COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-1

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITH S-9 ACTIVATION

<u>S. typhimurium</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: TA98 DATE PLATED: 08/18/2009 CELLS SEEDDED: 9.440E+07	REVERTANTS	973	23	18	16	18	17	26
	PER	758	28	22	21	22	23	19
	PLATE	934	23	17	22	17	16	24
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/18/2009 CELLS SEEDDED: 8.380E+07	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	862	76	50	97	83	120	100
	PER	932	78	55	90	86	96	92
	PLATE	1108	81	72	75	75	110	193
STRAIN: TA1535 DATE PLATED: 08/18/2009 CELLS SEEDDED: 1.824E+08	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	230	17	20	15	13	18	26
	PER	213	16	21	19	16	12	17
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDDED: 1.418E+08	PLATE	209	15	17	14	16	7	16
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	168	15	7	12	9	13	15
STRAIN: TA1537 DATE PLATED: 08/18/2009 CELLS SEEDDED: 1.418E+08	PER	157	10	13	8	23	14	15
	PLATE	172	10	10	8	9	13	14
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				100	500	1000	3000	5000
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDDED: 2.864E+08	REVERTANTS	488	12	20	14	16	16	21
	PER	527	18	18	13	22	33	31
	PLATE	478	20	16	20	16	20	25
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: WP2 uvrA DATE PLATED: 08/18/2009 CELLS SEEDDED: 2.864E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**CONFIRMATORY MUTATION ASSAY RAW COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-2

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

## WITHOUT ACTIVATION

<u><i>S. typhimurium</i></u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: TA98 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.078E+08	REVERTANTS	433	14	16	28	34	49	51
	PER	265	11	30	26	29	30	40
	PLATE	306	22	20	19	24	55	44
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.460E+08	REVERTANTS	256	59	49	62	96	114	260
	PER	229	87	42	65	81	93	196
	PLATE	218	79	45	49	76	129	141
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA1535 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.420E+08	REVERTANTS	525	18	12	15	15	29	26
	PER	469	15	16	11	9	13	14
	PLATE	515	17	18	18	19	21	23
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDED: 3.088E+08	REVERTANTS	178	4	4	5	6	5	6
	PER	225	11	5	2	6	6	8
	PLATE	193	11	4	6	4	4	3
	LAWN	NL	NL	NL	NL	NL	NL	NL
		PRECIPITATE	NP	NP	NP	NP	NP	NP

<u><i>E. coli</i></u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDED: 2.022E+08	REVERTANTS	409	8	15	15	10	17	17
	PER	473	12	9	19	13	17	8
	PLATE	440	14	23	19	17	17	14
	LAWN	NL	NL	NL	NL	NL	NL	NL
		PRECIPITATE	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY  
CONFIRMATORY MUTATION ASSAY RAW COLONY COUNTS AND BACKGROUND LAWN EVALUATION

EXPERIMENT NO.: B-2

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITH S-9 ACTIVATION

<u>S. typhimurium</u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: TA98 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.078E+08	REVERTANTS	959	25	30	24	30	32	55
	PER	956	31	29	21	25	33	46
	PLATE	844	40	21	24	31	32	41
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.460E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	972	106	82	68	92	102	161
	PER	984	82	130	72	80	118	193
	PLATE	774	79	92	60	103	96	326
STRAIN: TA1535 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.420E+08	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	152	13	28	15	29	10	16
	PER	177	11	23	15	21	9	12
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDDED: 3.088E+08	PLATE	225	12	19	26	19	18	11
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	225	8	10	5	8	10	3
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDDED: 3.088E+08	PER	258	8	6	8	13	5	5
	PLATE	218	6	4	13	2	9	7
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.022E+08	REVERTANTS	387	12	13	15	11	15	9
	PER	355	22	13	18	13	8	14
	PLATE	376	15	14	14	13	12	12
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.022E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**CONFIRMATORY MUTATION ASSAY CORRECTED COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-2

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITHOUT ACTIVATION

<u>S. typhimurium</u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: TA98 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.078E+08	REVERTANTS	461	18	20	33	39	55	57
	PER	283	15	35	31	34	35	45
	PLATE	327	26	24	23	28	61	50
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA100 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.460E+08	REVERTANTS	274	65	55	69	105	124	278
	PER	245	95	47	72	89	101	210
	PLATE	234	87	51	55	83	139	152
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1535 DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.420E+08	REVERTANTS	558	22	16	19	19	34	31
	PER	499	19	20	15	13	17	18
	PLATE	548	21	22	22	23	25	27
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDDED: 3.088E+08	REVERTANTS	191	7	7	8	9	8	9
	PER	241	15	8	5	9	9	12
	PLATE	207	15	7	9	7	7	6
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<u>E. coli</u>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDDED: 2.022E+08	REVERTANTS	435	12	19	19	14	21	21
	PER	503	16	13	23	17	21	12
	PLATE	468	18	27	23	21	21	18
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**SALMONELLA TYPHIMURIUM/ESCHERICHIA COLI PLATE INCORPORATION MUTATION ASSAY**  
**CONFIRMATORY MUTATION ASSAY CORRECTED COLONY COUNTS AND BACKGROUND LAWN EVALUATION**

EXPERIMENT NO.: B-2

SITEK STUDY NO.: 1003-2140

TEST ARTICLE: Ethylenediamine dinitrate (EDDN)

SOLVENT: DD Water

CONC. IN: µg/plate

WITH S-9 ACTIVATION

<i>S. typhimurium</i>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: TA98 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.078E+08	REVERTANTS	1017	30	35	28	35	37	61
	PER	1014	36	34	25	30	38	52
	PLATE	895	45	25	28	36	37	46
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: TA100 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.460E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	1031	115	90	75	100	111	173
	PER	1043	90	141	79	88	128	207
	PLATE	821	87	100	67	112	105	348
STRAIN: TA1535 DATE PLATED: 08/27/2009 CELLS SEEDED: 2.420E+08	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	164	17	33	19	34	14	20
	PER	190	15	27	19	25	13	16
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDED: 3.088E+08	PLATE	241	16	23	31	23	22	15
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	241	12	14	8	12	14	6
STRAIN: TA1537 DATE PLATED: 08/27/2009 CELLS SEEDED: 3.088E+08	PER	276	12	9	12	17	8	8
	PLATE	234	9	7	17	5	13	10
	LAWN	NL	NL	NL	NL	NL	NL	NL
	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP

<i>E. coli</i>		Positive Control	Solvent Control	Concentration per plate				
				3000	3500	4000	4500	5000
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDED: 2.022E+08	REVERTANTS	412	16	17	19	15	19	13
	PER	378	26	17	22	17	12	18
	PLATE	401	19	18	18	17	16	16
	LAWN	NL	NL	NL	NL	NL	NL	NL
STRAIN: WP2 uvrA DATE PLATED: 08/27/2009 CELLS SEEDED: 2.022E+08	PRECIPITATE	NP	NP	NP	NP	NP	NP	NP
	REVERTANTS	412	16	17	19	15	19	13
	PER	378	26	17	22	17	12	18
	PLATE	401	19	18	18	17	16	16

NL = Normal, healthy microcolony lawn.

NP = No precipitate.

**APPENDIX III**  
**SITEK'S HISTORICAL POSITIVE AND**  
**SOLVENT CONTROL DATA**

**SITEK RESEARCH LABORATORIES**

**HISTORICAL SOLVENT CONTROL DATA FOR SALMONELLA TYPHIMURUM/E. COLI  
PLATE INCORPORATION/PREINCUBATION MUTATION ASSAY  
MUTANT EXPRESSED PER PLATE  
WITHOUT S-4 ACTIVATION**

<b>TA99</b>	<b>DMSO</b>	<b>ACET</b>	<b>CORN OIL</b>	<b>H<sub>2</sub>O</b>	<b>SALINE</b>
AVERAGE	23	26	26	24	26
STANDARD DEVIATION (±)	5	5	3	4	6
MINIMUM VALUE	14	11	21	19	17
MAXIMUM VALUE	36	32	31	34	35
N*	52	20	12	34	25

<b>TA100</b>	<b>DMSO</b>	<b>ACET</b>	<b>CORN OIL</b>	<b>H<sub>2</sub>O</b>	<b>SALINE</b>
AVERAGE	57	99	77	64	74
STANDARD DEVIATION (±)	15	34	12	12	14
MINIMUM VALUE	36	26	51	45	50
MAXIMUM VALUE	132	174	118	99	132
N*	53	20	16	34	31

<b>TA1535</b>	<b>DMSO</b>	<b>ACET</b>	<b>CORN OIL</b>	<b>H<sub>2</sub>O</b>	<b>SALINE</b>
AVERAGE	13	19	15	14	15
STANDARD DEVIATION (±)	3	5	4	4	4
MINIMUM VALUE	9	12	11	8	10
MAXIMUM VALUE	20	33	23	29	24
N*	56	19	14	38	28

<b>TA1537</b>	<b>DMSO</b>	<b>ACET</b>	<b>CORN OIL</b>	<b>H<sub>2</sub>O</b>	<b>SALINE</b>
AVERAGE	9	9	10	8	9
STANDARD DEVIATION (±)	3	3	4	2	3
MINIMUM VALUE	2	6	3	5	4
MAXIMUM VALUE	20	16	23	12	15
N*	55	21	14	33	29

<b>E. COLI</b>	<b>DMSO</b>	<b>ACET</b>	<b>CORN OIL</b>	<b>H<sub>2</sub>O</b>	<b>SALINE</b>
AVERAGE	14	15	14	16	15
STANDARD DEVIATION (±)	3	4	3	3	4
MINIMUM VALUE	8	8	10	11	10
MAXIMUM VALUE	24	25	19	22	25
N*	53	19	13	35	28

N\* = NUMBER OF DATA POINTS.



# SITEK RESEARCH LABORATORIES

HISTORICAL SOLVENT CONTROL DATA FOR *SALMONELLA TYPHIMURUM* E. COLI  
 PLATE INCORPORATION/PREINCUBATION MUTATION ASSAY  
 MUTANT EXPRESSED PER PLATE  
 WITH S-9 ACTIVATION

<u>TA98</u>	<u>DMSO</u>	<u>ACET</u>	<u>CORN OIL</u>	<u>H<sub>2</sub>O</u>	<u>SALINE</u>
AVERAGE	80	31	31	33	32
STANDARD DEVIATION (±)	8	4	4	5	7
MINIMUM VALUE	10	23	22	22	20
MAXIMUM VALUE	48	37	35	43	52
N*	53	19	12	38	25

<u>TA100</u>	<u>DMSO</u>	<u>ACET</u>	<u>CORN OIL</u>	<u>H<sub>2</sub>O</u>	<u>SALINE</u>
AVERAGE	62	87	45	87	78
STANDARD DEVIATION (±)	18	30	12	14	15
MINIMUM VALUE	41	62	67	48	55
MAXIMUM VALUE	181	174	121	107	116
N*	53	19	16	37	31

<u>TA1535</u>	<u>DMSO</u>	<u>ACET</u>	<u>CORN OIL</u>	<u>H<sub>2</sub>O</u>	<u>SALINE</u>
AVERAGE	12	18	15	14	15
STANDARD DEVIATION (±)	3	6	3	3	3
MINIMUM VALUE	8	9	11	7	9
MAXIMUM VALUE	21	33	25	21	23
N*	56	19	14	38	23

<u>TA1537</u>	<u>DMSO</u>	<u>ACET</u>	<u>CORN OIL</u>	<u>H<sub>2</sub>O</u>	<u>SALINE</u>
AVERAGE	9	10	10	9	8
STANDARD DEVIATION (±)	3	3	4	3	2
MINIMUM VALUE	4	8	6	6	3
MAXIMUM VALUE	18	14	18	17	11
N*	57	19	22	34	29

<u>E. COLI</u>	<u>DMSO</u>	<u>ACET</u>	<u>CORN OIL</u>	<u>H<sub>2</sub>O</u>	<u>SALINE</u>
AVERAGE	16	17	18	17	19
STANDARD DEVIATION (±)	3	4	5	3	5
MINIMUM VALUE	10	12	0	10	10
MAXIMUM VALUE	23	26	26	24	27
N*	53	17	13	35	20

N\* = NUMBER OF DATA POINTS.

**SITEK RESEARCH LABORATORIES**

**HISTORICAL POSITIVE CONTROL DATA FOR *SALMONELLA TYPHIMURIUM*/E. COLI  
 PLATE INCORPORATION MUTATION/PREINCUBATION ASSAY  
 MUTANT EXPRESSED PER PLATE  
 WITH AND WITHOUT ACTIVATION**

<b>WITHOUT ACTIVATION</b>	<b>TA98 (2NF)</b>	<b>TA100 (NaAz)</b>	<b>TA1535 (NaAz)</b>	<b>TA1537 (9AA)</b>	<b>E. COLI (MMS)</b>
AVERAGE	618	419	338	135	444
STANDARD DEVIATION (±)	163	107	88	61	113
MINIMUM VALUE	180	233	57	24	122
MAXIMUM VALUE	1029	924	612	321	702
N*	107	107	109	107	108
<b>WITH ACTIVATION</b>	<b>TA98 (2AA)</b>	<b>TA100 (2AA)</b>	<b>TA1535 (2AA)</b>	<b>TA1537 (2AA)</b>	<b>E. COLI (2AA)</b>
AVERAGE	583	610	112	59	155
STANDARD DEVIATION (±)	324	261	49	28	83
MINIMUM VALUE	54	105	42	18	45
MAXIMUM VALUE	1732	1438	295	177	500
N*	311	110	109	110	108

N\* = NUMBER OF DATA POINTS.

**APPENDIX IV**

**STUDY PROTOCOL, PROTOCOL AMENDMENTS AND PROTOCOL DEVIATION**

**EVALUATION OF A TEST ARTICLE IN THE *SALMONELLA TYPHIMURIUM*/  
*ESCHERICHIA COLI* PLATE INCORPORATION MUTATION ASSAY IN THE  
PRESENCE AND ABSENCE OF INDUCED RAT LIVER S-9**

This protocol is presented in two parts. Part One is designed to collect specific information pertaining to the test article and study. Part Two describes the study design in detail. Please complete all bolded sections in Part One and sign section 8.0 to approve the protocol.

**PART ONE**

**1.0 SPONSOR**

1.1 Name: USA RDECOM, AMSRD-MSF

Environmental Acquisition & Logistics Sustaining Program

1.2 Address: Aberdeen Proving Ground, MD 21010

1.3 Sponsor's Study Coordinator: Gunda Reddy, Ph.D., DABT

**2.0 TESTING FACILITY**

2.1 Name: SITEK Research Laboratories

2.2 Address: 15235 Shady Grove Road, Suite 303  
Rockville, Maryland 20850

2.3 Study Director: Jian Song, Ph.D.

**3.0 STUDY NUMBERS**

\* 3.1 Testing Facility's Study No.: 1003-2140

3.2 Sponsor's Study No.: Not Available

**4.0 TEST ARTICLE**

GLP's require that test article characterization information must be provided in the final report. This includes identification, lot number, purity, stability, source, and expiration date. As per regulatory requirements, lack of the above information will be cited as a GLP violation in the "Study Director's Compliance Statement" section of the final report.

\* To be completed by the Testing Facility.

## 4.1 Identification

**Name:** Ethylenediamine dinitrate (EDDN)

**Batch/Lot No.:** To be provided

## 4.2 Description

**Color:** White to off-white

**Physical Form:** Powder

### 4.3 Analysis

**Purity Information:** 99.5%

**Does the Sponsor require the use of a correction factor to account for impurity?**

       **Yes**        X   **No**.

**If yes, what is the correction factor?** \_\_\_\_\_

**Determination of the test article characteristics as defined by Good Laboratory Practices will be the responsibility of the Sponsor. The specific GLP references for U.S. agencies are: FDA = 21 CFR, 58.105; EPA TSCA = 40 CFR, 792.105 and EPA FIFRA = 40 CFR 160.105.**

#### 4.4 Stability

**Storage Conditions (check one):**

  X   Dry/Room Temperature             Refrigerated (1-5°C)

### Frozen (-10 to -20°C)

**Other (please specify):** \_\_\_\_\_

**Expiration Date:** Not Available

**4.5 Preferred Solvent (check one):**

  X   H<sub>2</sub>O             DMSO             Acetone             Ethanol

**Other (please specify):** \_\_\_\_\_

**To be decided by the Testing Facility**

#### **4.6 Special Handling Instructions:**

**Use Standard Laboratory Safety Practices For Avoiding Exposure To  
Hazardous Substances And Follow Safety Requirement For Explosive Material.**

### **5.0 REGULATORY AGENCY SUBMISSION**

#### **5.1 Test Design Specifications**

This study protocol is designed to meet or exceed the US EPA, ICH and OECD Guidelines specified in the following documents (1, 2, 3):

United States Environmental Protection Agency, Title 40 Code of Federal Regulations, Part 798, Health Effects Testing Guidelines, Subpart F, Sec. 798.5265, the *Salmonella typhimurium* reverse mutation assay. Revised July 1, 2002.

OECD Guideline for the Testing of Chemicals, No. 471. Bacterial Reverse Mutation Test. Revised July 21, 1997.

International Conference on Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use. ICH Harmonised Tripartite Guideline S2A. Guidance on Specific Aspects of Regulatory Genotoxicity Tests for Pharmaceuticals. Federal Register 61 (80): 18198-18202, 1996.

#### **5.2 Good Laboratory Practices**

This study will be conducted in compliance with the following Good Laboratory Practice standards:

United States Environmental Protection Agency, Title 40 Code of Federal Regulations Parts 160 and 792. Revised July 1, 2002.

United States Food and Drug Administration, Title 21 Code of Federal Regulations Part 58. Revised April 1, 2003.

Japanese Ministry of Agriculture, Forestry and Fisheries, 11 NohSan, Notification No. 6283, October 1, 1999.

Japanese Ministry of Health and Welfare, Ordinance No. 21, April 1, 1997.

Japanese Ministry of International Trade and Industry, Notification No. 85, Basic Industries Bureau, March 31, 1984.

Organisation for Economic Cooperation and Development, The OECD Principles of Good Laboratory Practice, Environment Monograph No. 45 [ENV/MC/CHEM(98)17], Paris 1998.

Will this study be submitted to a regulatory agency?

☒ Yes      ☐ No

If so, which agency(ies)? Worldwide

## 6.0 TEST ARTICLE/DOSING SOLUTIONS CHARACTERIZATION

The U.S. requirements for analysis of dosing solutions are specified in: FDA = 21 CFR, 58.113; EPA TSCA = 40 CFR, 792.113; and EPA FIFRA = 40 CFR, 160.113.

Does the Sponsor want dosing solution analysis?

☐ Yes\*\*      ☒ No

If yes, please complete the rest of this section.

If requested by the Sponsor, SITEK Research Laboratories will determine the strength and stability of the dosing solutions. The method of analysis may be provided by the Sponsor, or if requested by the Sponsor, SITEK Research Laboratories will develop the method of analysis.

Alternatively, the Sponsor will be responsible for determining the strength and stability of the dosing solutions.

Dosing solution analysis will be performed by:

☐ SITEK Research Laboratories      ☐ Sponsor\*\*\*

What dosing solutions will be analyzed? \_\_\_\_\_

---

\*\* Additional charges will apply. See Special Services price schedule.

\*\*\* Please note: All work pertaining to this study that is performed outside of SITEK is the responsibility of SITEK's Study Director. Therefore, as required by the GLPs, all of the following must be forwarded to the Study Director:

- All subcontract and/or Sponsor Quality Assurance audit findings and comments.
- Any deviations and/or amendments, if applicable.
- An original or copy of the analysis report.
- Location of where the raw data from the analysis will be archived.

If the subcontract work is not performed under the GLPs, a statement by the Sponsor informing SITEK's Study Director of such must be provided.

**From the Range Finding Test?**

\_\_\_\_\_ Yes                      \_\_\_\_\_ No

**From the Assay?**

\_\_\_\_\_ Yes                      \_\_\_\_\_ No

**Which concentration(s)?** \_\_\_\_\_

**What amount of each concentration?** \_\_\_\_\_

**At what temperature should the dosing solutions be stored?**

\_\_\_\_\_ Room Temperature                      \_\_\_\_\_ Frozen (-10 to -20° C)  
\_\_\_\_\_ Refrigerated (1-5° C)

**At what temperature should the dosing solutions be shipped?**

\_\_\_\_\_ Room Temperature                      \_\_\_\_\_ On Wet Ice  
\_\_\_\_\_ On Dry Ice

**7.0 STUDY DATES**

\* 7.1 Proposed Experimental Start Date: August 7, 2009

Defined as the first date the test article is applied to the test system.

\* 7.2 Anticipated Experimental Completion Date: September 18, 2009


Defined as the last date on which data are collected directly from the study.

\* 7.3 Anticipated Draft Report Submission Date: September 25, 2009

7.4 Final Report: The final report will be initiated sixty days after remittance of the draft report and issued no later than thirty days thereafter.



**8.0 PROTOCOL APPROVAL**

\*   
\_\_\_\_\_  
Study Director


8-5-2009  
Date

  
\_\_\_\_\_  
Sponsor's Study Coordinator

8-5-2009  
Date

\*   
\_\_\_\_\_  
Quality Assurance Manager

8-5-09  
Date

\*   
\_\_\_\_\_  
Safety Officer

8-5-09  
Date

\* To be completed by the Testing Facility.

## STUDY DESIGN

### PART TWO

#### 9.0 PURPOSE

The purpose of this study is to evaluate the test article for its potential to cause mutations in the histidine operon of *Salmonella typhimurium* strains TA98, TA100, TA1535 and TA1537 and the tryptophan operon of *Escherichia coli* strain WP2 uvrA.

#### 10.0 JUSTIFICATION FOR SELECTION OF TEST SYSTEM

The *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay has been used extensively and has been demonstrated to be effective in detecting the mutagenic activity of chemicals from a wide range of classes.

#### 11.0 ABBREVIATIONS

2-AA	-	2-Aminoanthracene
2-NF	-	2-Nitrofluorene
9-AA	-	9-Aminoacridine
DMSO	-	Dimethyl Sulfoxide
MMS	-	Methyl Methanesulfonate
NaN <sub>3</sub>	-	Sodium Azide
NADP	-	Nicotinamide-adenine Dinucleotide Phosphate
O.D.	-	Optical Density
%T	-	Percent Transmittance
S-9	-	Induced Rat Liver Homogenate

#### 12.0 INDICATOR CELLS

##### 12.1 Source

The *Salmonella typhimurium* strains TA98, TA100, TA1535 and TA1537 were obtained from Dr. Bruce N. Ames, University of California, Berkeley, California. The *Escherichia coli* strain WP2 uvrA was obtained from Ms. Judy Mayo of Pharmacia Corporation, Kalamazoo, Michigan.

## 12.2 Culture Conditions

The *Salmonella typhimurium* and *Escherichia coli* strains are routinely grown in Oxoid Nutrient Broth No. 2 in a shaker incubator rotating at approximately 120 rpm and maintaining a temperature of  $37 \pm 1^\circ\text{C}$ .

## 12.3 Stock Cultures

The *Salmonella typhimurium* and *Escherichia coli* strains were propagated to obtain a sufficient number of cells for freezing a large number of stock ampules. The cells were cryopreserved in Oxoid Nutrient Broth No. 2 supplemented with 8-9% dimethyl sulfoxide (DMSO) and stored in liquid nitrogen vapor phase. Scrapes from stock ampules are used to initiate the stock cultures for the test.

## 13.0 METABOLIC ACTIVATION

The standard rat liver S-9 prepared from male Sprague-Dawley rats with Aroclor-1254 or Phenobarbital and/or  $\beta$ -naphthoflavone will be used for the metabolic activation system

## 14.0 ROUTE OF ADMINISTRATION OF TEST ARTICLE

The test article will be administered in vitro directly or through a solvent compatible with the test cultures. This is the only route of administration available in this test system.

## 15.0 TEST SYSTEM IDENTIFICATION

All test plates will be labeled using an indelible pen with a code system which clearly identifies the experiment number, the SITEK test article number, controls, doses, and whether or not the plate was treated in conjunction with an exogenous activation system.

The test article will be designated by the unique four-digit number assigned by SITEK when the test article is received (e.g., 0074). The experiment phase will be designated by the letter A (Range Finding Test) or B (Mutation Assay) followed by a number designating the trial number. This will be followed by the letter N (No Activation) or S (With S-9) which will be followed by the dose and strain identification numbers. The doses will be identified by the numbers 1, 2, 3, ... indicating the highest to the lowest dose. The strain identification numbers will be as follows:

*Salmonella typhimurium*      *Escherichia coli*

1 = TA98	5 = WP2 uvrA
2 = TA100	
3 = TA1535	
4 = TA1537	

An example of a plate label follows:

0074B1-S-1-3

0074	=	SITEK Test Article Number
B1	=	First Mutation Assay
S	=	With S-9
1	=	Highest Test Article Dose
3	=	Strain TA1535

In addition to the above, the Range Finding Test and Mutation Assay viability plates that contain 10X (0.5mM) histidine biotin or 10X (0.5mM) tryptophan will be designated with the prefix "T".

## 16.0 CONTROL SUBSTANCES

### 16.1 Positive Controls

The positive control chemicals that will be used for the tester strains in the presence and absence of exogenous metabolic activation are presented below. The abbreviations are defined in Section 11.0.

	<u>Strain</u>	<u>S-9</u>	<u>Chemical</u>	<u>Dose</u> <u>(µg/plate)</u>
<i>Salmonella typhimurium</i>				
	TA98	-	2-NF	2.5-7.5
	TA98	+	2-AA	1.25-5.0
	TA100	-	NaAz	0.5-2.0
	TA100	+	2-AA	1.25-5.0
	TA1535	-	NaAz	0.5-2.0
	TA1535	+	2-AA	1.25-5.0
	TA1537	-	9-AA	25-75
	TA1537	+	2-AA	1.25-5.0

### *Escherichia coli*

WP2 <i>uvrA</i>	-	MMS	2000-4000
WP2 <i>uvrA</i>	+	2-AA	10-20

If necessary, other appropriate positive controls can be used with the approval of the Sponsor.

DMSO will be used to solubilize the positive controls, except for NaAz and MMS, which will be dissolved in deionized, distilled H<sub>2</sub>O.

#### 16.2 Solvent Control

The solvent used for dissolving the test article will be used as the solvent control. Deionized, distilled water, dimethyl sulfoxide (CAS #67-68-5), ethanol (CAS #64-17-5) and acetone (CAS #67-64-1) are some of the solvents which are compatible with this test system. If there is a need to use other solvents, the approval of the Sponsor will be obtained prior to their use.

### 17.0 DOCUMENTATION

All procedures, results, significant observations, and methods used for analysis of results will be documented in a study notebook. The study notebook will also include copies of the protocol, all protocol amendments and protocol deviations, study reports, and all relevant communications with the Sponsor.

### 18.0 EXPERIMENTAL PROCEDURE

#### 18.1 Determination of Solubility/Miscibility

In order to determine the optimal vehicle for delivering the test article to the test system or to determine the maximum achievable concentration in the solvent requested by the Sponsor, a solubility/miscibility test may be performed, if necessary. The solvents of choice for this system are water, DMSO, acetone and ethanol. If the test article is not sufficiently soluble in any of these solvents, additional solvents will be screened.

For solid and viscous test articles, the solubility test will consist of weighing out 20- to 100-mg aliquots of test article and adding solvent in 0.1 mL increments, with thorough mixing between additions, until the test article is dissolved as determined by visual inspection or until 5.0 mL of solvent has been added to the vessel. The volume of solvent required for complete dissolution and any additional observations will be recorded in the study notebook. Test articles that do not dissolve in 5.0 mL of solvent will be visually inspected and recorded as either "not soluble," "partially soluble forming a homogeneous suspension," or "partially soluble not forming a homogeneous suspension."

For liquid test articles, a miscibility test will be conducted. 0.5 mL of solvent will be added to 0.5 mL aliquots of the test article. The resulting solution will be thoroughly mixed and observed for miscibility. The test article will be rated by visual inspection as either "not miscible," "partially miscible," or "completely miscible" in each of the four preferred solvents. The miscibility rating and any additional observations will be recorded in the study notebook.

Where solubility cannot be achieved, the test article will be delivered as a suspension in the desired vehicle. If sufficient solubility data is available, the solubility/miscibility test will not be performed.

## 18.2 Preparation of Test Cultures

The strains of *Salmonella typhimurium* and *Escherichia coli* will be prepared from cultures that were started from scrapes placed in Oxoid Nutrient Broth No. 2. The cultures will be placed on the shaker, and a timer turns on the incubator approximately 8-12 or 4-6 hours for *Salmonella typhimurium* or *Escherichia coli*, respectively, prior to sampling the cultures for growth determination. The incubator will be set at 120 rpm and  $37 \pm 1^\circ\text{C}$ . Samples from each culture will be checked for Percent Transmittance (%T) at 650 nm.

Only cultures that have a %T of between 25% (O.D. 0.6) and 10% (O.D. 1.0) will be used.

## 18.3 Preparation of S-9 Metabolic Activation Mix

For the portion of the Range Finding Test or the Mutation Assay in which the cells are exposed to the test article in conjunction with an exogenous metabolic activation system, induced rat liver S-9 plus cofactors (S-9 mix) will be used as the activation system. The components of the standard S-9 mix will be 8mM  $\text{MgCl}_2$ , 33mM KCl, 5mM glucose-6-phosphate, 4mM NADP, 100mM sodium phosphate buffer (pH 7.4), and 10% rat liver S-9.

## 18.4 Preparation of Test Article

The desired amount of the test article as specified in the dilution scheme will be weighed or measured just prior to use in either the Range Finding Test or the Mutation Assay. The dosing solutions will be prepared by adding the appropriate volume of solvent to the test article and thoroughly mixing the resulting solution until the test article goes completely into solution or a homogeneous suspension is achieved. The remaining doses specified in the dilution scheme will be prepared by either performing a serial dilution or by varying the volume delivered from the stock concentration to the cultures. In all treatments the amount of solvent delivered to the target cultures will be limited to a level which has no cytotoxic effect on the cells. If necessary, the test article may be added directly to the top agar.

## 18.5 Range Finding Test

In order to determine the test article concentrations that will produce from 0-100% toxicity, a Range Finding Test will be performed with and without S-9 activation using tester strains TA100 and WP2 uvrA only. The test article will be weighed or measured, and a serial dilution will be prepared. If there are no solubility/miscibility limitations, prior knowledge of cytotoxicity indicates differently, or the Sponsor specifies differently, the treatment concentrations for solid and viscous test articles will be 5000, 1000, 500, 100, 50, 10 and 5.0  $\mu\text{g}/\text{plate}$ . If the results based on the dosing regimen indicate that the threshold level of complete toxicity is below 5.0  $\mu\text{g}/\text{plate}$  an additional Range Finding Test will be performed.

### 18.5.1 Treatment

2.0 mL aliquots of molten top agar, to which trace amounts of histidine and biotin have been added, will be dispensed to a series of culture tubes maintained at  $45 \pm 1^\circ\text{C}$ . Treatment will be performed by adding 0.5 mL of S-9 mix or 0.5 mL of sterile, distilled, deionized water, 0.1 mL of tester strain TA100 or WP2 uvrA, and 0.1 mL of test article to the top agar. Appropriate solvent controls will also be prepared.

In addition, plates for determining viability will be prepared by plating the test article doses with a  $2.0 \times 10^5$  dilution of tester strain TA100 or WP2 uvrA in top agar containing 10X histidine-biotin or 10X tryptophan, respectively.

The contents will be mixed by vortexing the tube, and then the contents will be poured onto a bottom agar plate and evenly distributed by gently tilting and rotating the plate. The plate will be placed on a flat, level surface until solidified. After all treatment is performed, the plates will be inverted and incubated at  $37 \pm 1^\circ\text{C}$  for 48-72 hours.

#### 18.5.2 Determination of Toxicity

After 48-72 hours of incubation, the plates will be removed from the incubator and evaluated or placed in cold storage ( $1-5^\circ\text{C}$ ) until evaluated.

Evaluation of test article toxicity on the tester strain will be based on three end points:

1. Viability of cells plated on minimal medium plates supplemented with excess histidine-biotin or tryptophan. Toxicity will be measured as a decrease in the number of colonies per plate with increasing test article concentration.
2. The number of revertant colonies on minimal medium plates supplemented with trace amounts of histidine-biotin or tryptophan. Toxicity will be measured as a reduction in the number of revertant colonies per plate with increasing test article concentration.
3. The integrity of the background microcolony lawn. Toxicity will be measured as a thinning or disappearance of the background lawn usually occurring with an increase in the size of the remaining microcolonies relative to the control plates.

The number of revertants per plate and the number of viable colonies per plate will be determined by counting them with an automatic colony counter or by hand as described in Sections 18.6.5.1 and 18.6.5.2.

The counts will be entered directly in the Excel 2003 computer program 2140A, and the calculations will be performed. The computer printouts will be included in the study notebook.

#### 18.6 Mutation Assay

The maximum concentration of nontoxic test articles that is tested will be 5 mg per plate, unless the Sponsor requests otherwise or precipitation of the test article on the plate warrants the use of a lower concentration. Test articles that produce a toxic effect will be tested at a maximum dose that significantly reduces the number of revertants per plate and/or causes thinning of the background lawn. Four lower doses will be selected that should not produce toxicity. Test articles that are insoluble at concentrations of 5 mg per plate or lower will be tested at a maximum dose that produces precipitate. A concentration that produces precipitate in the test system will be considered to be beyond the limits of solubility. The actual dose levels for the assay, once determined, will be added to the protocol in the form of an amendment. Each test article dose, the positive controls and solvent controls will be plated in triplicate.

#### 18.6.1 Test Culture Preparation and Exposure

Cultures of *Salmonella typhimurium*, TA98, TA100, TA1535, TA1537, and *Escherichia coli* WP2 uvrA for use in the Mutation Assay will be prepared as described in Section 18.2. The test article will be weighed or measured, and a serial dilution will be performed as previously described in Section 18.4. 2 mL aliquots of molten top agar to which histidine and biotin or tryptophan have been added will be dispensed to a series of culture tubes maintained at  $45 \pm 1^\circ\text{C}$ . Treatment will be performed by adding 0.5 mL of S-9 mix or 0.5 mL of sterile, distilled, deionized water, 0.1 mL of tester strain, and 0.1 mL of test article to the top agar. Appropriate solvent and positive controls will also be prepared. The contents will be mixed by vortexing the tube, and then the contents will be poured onto a bottom agar plate and evenly distributed by gently tilting and rotating the plate. The plate will be placed on a flat, level surface until solidified. After all treatment will be performed, the plates will be inverted and incubated at  $37 \pm 1^\circ\text{C}$  for 48-72 hours.

#### 18.6.2 Confirmation of Tester Strain Genotypes

On the same day as the plating of the Mutation Assay, the genotypes of the tester strains will be confirmed. All of the *Salmonella typhimurium* strains will be tested for histidine dependence and the rfa mutation. Each *Salmonella typhimurium* strain will be tested for the uvrB deletion after cryopreservation of the stock ampules. The tester strains TA98 and TA100 will also be tested for the pKM101 plasmid. The *Escherichia coli* WP2 uvrA strain will be tested for tryptophan dependence.

#### 18.6.3 Tester Strain Viability Determination

After the Mutation Assay has been plated, a dilution of each tester strain will be prepared, and approximately 250-500 bacteria will be plated in top agar supplemented with 10X histidine-biotin or 10X tryptophan. These plates will be incubated for 48-72 hours, and then the total number of colonies that develop will be determined.

#### 18.6.4 Background Lawn Evaluation

The integrity of the background microcolony lawn will be evaluated by viewing each plate with the aid of a 2X to 4X microscope. The lawns will be rated as normal, slightly reduced, markedly reduced, extremely reduced or absent.

#### 18.6.5 Enumeration of Colonies

After 48-72 hours of incubation, the plates treated with the highest test article concentration will be observed for the presence of precipitate. If precipitate is absent, the entire assay will be counted using an automatic colony counter. If observation of the high dose plates reveals precipitate that interferes with accurate automatic counting, those plates will be counted by hand. The procedure will be repeated for each subsequent dose level or until no precipitate is evident.



#### 18.6.5.1 Automatic Colony Counting

Each plate will be placed on the stage, and three counts are made with the automatic counter. The plate will be rotated on the stage approximately 120° between each count, and the median count will be recorded.

#### 18.6.5.2 Hand Counting

Hand counting of colonies will be performed by marking a dot over each colony on the bottom of the plate while clicking off the counts on a digitometer. The hand count will be recorded for each plate.

The counts will be entered directly in the Excel 2003 computer program 2140B. The computer printouts will be included in the study notebook.

#### 18.7 Confirmatory Mutation Assay

If the first Mutation Assay gives negative or equivocal results, a confirmatory Mutation Assay will be performed. The test article treatment concentrations may be altered based on the results obtained in the first Mutation Assay. On the other hand, if the results of the first Mutation Assay are clearly positive, a confirmatory Mutation Assay may or may not be performed depending on the Sponsor's instructions.

#### 18.8 Criteria For a Valid Assay

The following criteria will be used as guidelines in determining the acceptability of the results. Since it is impossible to formulate criteria that would apply to every configuration of data generated by the Mutation Assay, the Study Director will be responsible for the ultimate decision regarding the acceptability of the results.

##### 18.8.1 Solvent Control Cultures

The mean reversion frequency of the test article solvent control plates for each strain must fall within the range presented below.

*Salmonella typhimurium* *Escherichia coli*

TA98	30 ± 15	WP2 uvrA	15 ± 10
TA100	100 ± 70		
TA1535	20 ± 15		
TA1537	15 ± 12		

##### 18.8.2 Positive Controls

The results for the positive control cultures will be considered acceptable if the treated strains have mean reversion frequencies that are three times or greater than the mean reversion frequencies of the test article solvent control plates.

### 18.8.3 Tester Strain Characterization

1. All of the *Salmonella typhimurium* strains will be confirmed positive for histidine dependence and the *Escherichia coli* strain for tryptophan dependence.
2. All of the *Salmonella typhimurium* strains will be confirmed positive for the *rfa* mutation as evidenced by sensitivity to crystal violet.
3. The R-factor strains, TA98 and TA100, will be confirmed positive for the pKM101 plasmid as evidenced by ampicillin resistance.
4. The titer of the stock cultures of each strain will indicate that the stock cultures contained greater than  $0.5 \times 10^9$  cells/mL.

### 18.9 Evaluation of Test Results

The following criteria will be used as guidelines in evaluating the results of the Mutation Assay for a negative, positive or equivocal response. Since it is impossible to write criteria that would apply to every configuration of data generated by the Mutation Assay, the Study Director will be responsible for the ultimate decision in the evaluation of the results. The factors considered in making the decision will be discussed in the report.

#### 18.9.1 Criteria for a Negative Response

A response will be considered negative if 1) strains TA98 and TA100 have mean reversion frequencies that are less than twice that of the mean reversion frequencies of the corresponding solvent control plates, 2) strains TA1535, TA1537 and WP2 *uvrA* have mean reversion frequencies less than three times that of the corresponding solvent control plates, and 3) there is no evidence of a dose-dependent response.

#### 18.9.2 Criteria for a Positive Response

A response will be considered positive if either strain TA98 or TA100 has a dose that produces a mean reversion frequency that is greater than or equal to two times the mean reversion frequency of the corresponding solvent control plates or if either strain TA1535, TA1537 or WP2 *uvrA* has a dose producing a three-fold or greater increase in the mean reversion frequency compared to the solvent control frequency. In addition, the response must be dose-dependent or increasing concentrations of the test article must show increasing mean reversion frequencies. In evaluating the results, consideration will be given to the degree of toxicity exhibited by the dose causing the two-fold/three-fold or greater increase in reversion frequency and the magnitude of the increase in reversion frequency.

#### 18.9.3 Criteria for an Equivocal Response

A response will be considered equivocal if it does not fulfill the criteria of either a negative or a positive response and/or the Study Director does not consider the response to be either positive or negative.

In addition, if either strain TA1535, TA1537 or WP2 uvrA has a dose producing a twofold increase in mean reversion frequency compared to the solvent control frequency and there is a dose-dependent response at lower concentrations in this strain, then this result will be considered equivocal and the test may be repeated after consultation with the Sponsor.

## 19.0 PROTOCOL AMENDMENTS AND DEVIATIONS

If changes in the approved protocol are necessary, such changes will be documented in the form of protocol amendments and protocol deviations. Protocol amendments will be generated when changes in the protocol are made prior to performing a study or part of a study affected by the changes. In such cases, a verbal agreement to make such changes will be made between the Study Director and the Sponsor. These changes and the reasons for them will be documented and attached to the protocol as an addendum. Protocol deviations will be generated when the procedures used to perform the study do not conform to the approved protocol. The Sponsor will be informed of these deviations, and as soon as practical, such changes along with their reasons or explanations will be documented and kept in the study notebook.

## 20.0 REPORT OF RESULTS

### 20.1 Content

The results of the study will be submitted to the Sponsor in the form of a final report. A draft report will be submitted before the final report is issued. The final report will be initiated sixty days after remittance of the draft report and issued no later than thirty days thereafter. The report will include, but not be limited to, the following:

1. Name and address of the testing facility and the dates on which the study was initiated and completed, terminated or discontinued.
2. Objectives and procedures stated in the approved protocol, including any changes in the original protocol.
3. Methods used to analyze the data.
4. The test and control substances.
5. Description of the methods used to perform the study.
6. The data, mean plate counts, +/- SD, and any observations regarding toxicity and precipitate.
7. The name and signature of the Study Director and the names of other technical personnel who participated in performing the study.
8. The location where the raw data and reports are to be stored.
9. A statement from the Quality Assurance Unit.

## 20.2 Changes and Corrections to the Final Report

All changes to the final report will be in the form of a report amendment which will include the reason(s) for the change, and the amendment will be added to the final report as an addendum.

## 21.0 ARCHIVES

The raw data, electronic file containing the data tables, documentation, protocol and final report of the study will be maintained in the SITEK Research Laboratories Archives, 15235 Shady Grove Road, Suite 303, Rockville, Maryland, according to the terms and conditions of the study.

## 22.0 REFERENCES

1. Ames, B. N., J. McCann and E. Yamasaki. Methods for detecting carcinogens and mutagens with the *Salmonella*/mammalian-microsome mutagenicity test. *Mut. Res.*, 31:347-367, 1975.
2. Maron, D., and B. N. Ames. Revised methods for the *Salmonella* mutagenicity test. *Mut. Res.*, 113:173-215, 1983.
3. Green, M. H. L., and W. J. Muriel. Mutagen testing using trp<sup>+</sup> reversion in *Escherichia coli*. In: B. J. Kilbey, et al. (eds.), *Handbook of Mutagenicity Test Procedures*, pp. 65-94, Elsevier North Holland Biomedical Press, Amsterdam, 1977.
4. Venitt, S., and J. M. Parry (eds.). *Mutagenicity testing: A practical approach*. IRL Press, Oxford, England and Washington, D.C., 1984.

**PROTOCOL AMENDMENT**

Amendment No.: 1

Sponsor: USA RDECOM, AMSRD-MSF,  
Environmental Acquisition & Logistics Sustaining  
Program,  
Aberdeen Proving Ground, MD 21010

Testing Facility: SITEK Research Laboratories  
15235 Shady Grove Road, Suite 303  
Rockville, Maryland 20850

SITEK's Study No.: 1003-2140

Sponsor's Study No.: N/A

Test Article I.D.: Ethylenediamine dinitrate (EDDN)

Protocol Title: Evaluation of a Test Article in the *Salmonella*  
*Typhimurium*/*Escherichia Coli* Plate Incorporation  
Mutation Assay in the Presence and Absence of  
Induced Rat Liver S-9

Amendment No. 1: Protocol page 12, section 18.6: Based on the results of  
the Range Finding Test the actual dose levels for the  
Definitive Mutation Assay were 100, 500, 1000, 3000  
and 5000 µg/plate for *Salmonella Typhimurium* and  
*Escherichia Coli* both with and without activation. The  
actual dose levels for the Confirmation Mutation Assay  
were 3000, 3500, 4000, 4500 and 5000 µg/plate for  
*Salmonella Typhimurium* and *Escherichia Coli* both with  
and without activation.

Reason for Amendment No. 1: Protocol page 12, section 18.6: The actual dose levels for  
the assay, once determined, will be added to the protocol  
in the form of an amendment.

APPROVED:



Jian Song, Ph.D.  
Study Director

8-21-09

Date

**PROTOCOL AMENDMENT**

Amendment No.: 2

Sponsor: USA RDECOM, AMSRD-MSF  
Environmental Acquisition & Logistics  
Sustaining Program  
Aberdeen Proving Ground, MD 21010

Testing Facility: SITEK Research Laboratories  
15235 Shady Grove Road, Suite 303  
Rockville, Maryland 20850

SITEK's Study No.: 1003-2140

Sponsor's Study No.: N/A


Test Article ID: Ethylenediamine dinitrate (EDDN)

Protocol Title: Evaluation of a Test Article in the *Salmonella typhimurium*/*Escherichia coli* Plate Incorporation Mutation Assay in the Presence and Absence of Induced Rat Liver S-9

Amendment No. 2: Protocol Page 1, Section 2.3, Study Director, Jian Song, Ph.D. has been replaced by Paul E. Kirby, Ph.D. as Study Director.

Reason for Amendment No. 2: Jian Song, Ph.D. is no longer in the employ of SITEK Research Laboratories.

**APPROVED:**

  
\_\_\_\_\_  
Paul E. Kirby, Ph.D.  
Study Director

2-22-10  
\_\_\_\_\_  
Date

**PROTOCOL DEVIATION**

Deviation No.: 1

Sponsor: USA RDECOM, AMSRD-MSF,  
Environmental Acquisition & Logistics Sustaining  
Program,  
Aberdeen Proving Ground, MD 21010

Testing Facility: SITEK Research Laboratories  
15235 Shady Grove Road, Suite 303  
Rockville, Maryland 20850

SITEK's Study No.: 1003-2140

Sponsor's Study No.: N/A

Test Article I.D.: Ethylenediamine dinitrate (EDDN)

Protocol Title: Evaluation of a Test Article in the *Salmonella*  
*Typhimurium*/*Escherichia Coli* Plate Incorporation  
Mutation Assay in the Presence and Absence of  
Induced Rat Liver S-9.

Deviation No. 1: Protocol page 13, section 18.6.1: Test Culture  
Preparation and Exposure: After all treatment will be  
performed, the plates will be inverted and incubated at  $37 \pm 1$  °C for 48-72 hours. However, the warm room  
temperature was from 34 to 36 °C during the study.

Reason for Deviation No. 1: The minimum warm room temperature was 34 °C  
because the door of the warm room was opened several  
times during the study and the occasion low temperature  
does not affect the result of the assays.

APPROVED:



Jian Song, Ph.D.  
Study Director

9-28-09

Date

**APPENDIX V**

**S-9 BATCH INFORMATION**



## MOLTOX POST MITOCHONDRIAL SUPERNATANT (S-9) QUALITY CONTROL & PRODUCTION CERTIFICATE

LOT NO.: <u>2342</u>	SPECIES: <u>Rat</u>	PREPARATION DATE: <u>September 9, 2008</u>
PART NO.: <u>11-101</u>	STRAIN: <u>Sprague Dawley</u>	EXPIRATION DATE: <u>September 9, 2010</u>
VOLUME: <u>2 ml</u>	SEX: <u>Male</u>	BUFFER: <u>0.154 M KCl</u>
	TISSUE: <u>Liver</u>	INDUCING AGENT(s): <u>Aroclor 1254</u>
REFERENCE: <u>Maron, D &amp; Ames, B, <i>Mutat Res</i> 113:173, 1983</u>		<u>(Monsanto KL615), 500 mg/kg i.p.</u>
STORAGE: <u>At or below -70°C</u>		

**BIOCHEMISTRY:****- PROTEIN**35.5 mg/mlAssayed according to the method of Lowry et al., *JBC* 193:265, 1951 using bovine serum albumin as the standard.**- ALKOXYRESORUFIN-0-DEALKYLASE ACTIVITIES**

Activity	P450	Fold - Induction
EROD	1A1, 1A2	72.8
PROD	2B1, 2B2	14.3
BROD	2B1, 2B2	27.2
MROD	1A1, 1A2	30.4

Assays for ethoxyresorufin-0-deethylase (EROD), pentoxy-, benzyl- and methoxyresorufin-0-dealkylases (PROD, BROD, & MROD) were conducted using a modification of the methods of Burke, et al., *Biochem Pharm* 34:3337, 1985. Fold-inductions were calculated as the ratio of the sample vs. uninduced specific activities (SA's). Control SA's (pmoles/min/mg protein) were 39.0, 17.4, 49.8, & 16.3 for EROD, PROD, BROD and MROD, respectively.

**BIOASSAY:****- TEST FOR THE PRESENCE OF ADVENTITIOUS AGENTS**

Samples of S-9 were assayed for the presence of contaminating microflora by plating 1.0 ml volumes on Nutrient Agar and Minimal Glucose (Vogel-Bonner E, supplemented with 0.05 mM L-histidine and D-biotin) media. Triplicate plates were read after 40 - 48 h incubation at 35 ± 2°C. The tested samples met acceptance criteria.

**- PROMUTAGEN ACTIVATION**

No. His+ Revertants

TA98    TA1535

171.6    960

The ability of the sample to activate ethidium (EtBr) EtBr/CPA and cyclophosphamide (CPA) to intermediates mutagenic to TA98 and TA1535, respectively, was determined according to Lesca, et al., *Mutation Res* 129:299, 1984. Data were expressed as revertants per µg EtBr or per mg CPA.

Dilutions of the sample S9, ranging from 0.2 - 10% in S9 mix, were tested for their ability to activate benzo(a)pyrene (BP) and 2-aminoanthracene (2-AA) to intermediates mutagenic to TA100. Assays were conducted as described by Maron & Ames, (*Mutat Res* 113:173, 1983).

µl S9 per plate/number his<sup>+</sup> revertants per plate

Promutagen	0	1	5	10	20	50
BP (5 µg)	152	286	295	499	653	993
2-AA (2.5 µg)	155	288	915	1563	2376	2064


**MOLECULAR TOXICOLOGY, INC.**

157 Industrial Park Dr.

Boone, NC 28607

(828) 264-9099

www.moltox.com




## CERTIFICATE OF ANALYSIS

Product Name: **Rat liver Postmitochondrial Supernatant (S-9)**

Product Number: 74002  
 Lot Number: 1102-2001  
 Product Brand: APSciences, Inc.  
 Storage Temperature: -80°C

Strain	Sprague Dawley	
Species	Rat	
Tissue Type	Liver	
Appearance	Reddish brown liquid in 2mL vial	
Expiration Date	2 Years	June 19, 2011
Protein Concentration	38.5 mg/mL in 0.15 KCl	
Preparation	Inducing agents: Phenobarbital (75 mg/kg body weight); b-naphthoflavone (80 mg/kg body weight); intraperitoneal injection once per day from days 1 to 4 (phenobarbital), and from day 3 to day 4 (b-naphthoflavone). Organs harvested for S9 preparation on day 5.	
BioAssays	<ol style="list-style-type: none"> <li>1. Chromosome Aberration Assay in Human Peripheral Blood Lymphocytes: Induction of cells with chromosomal aberrations observed with cyclophosphamide (no aberrant cells in solvent control).</li> <li>2. Salmonella typhimurium and E. Coli WP2 uvrA Mutation Assays: Induction of revertants by 2-aminoanthracene (2.5 ug/plate) in TA98 (50 x background); TA100 (8 x background); TA1535 (11 x background); TA1537 (9 x background); and 2-amininoanthracene (20 ug/plate) in WP2 uvrA (6 x background).</li> </ol>	

Signature:   
 Printed Name: Nola Mahaney  
 Title: VP, Operations  
 Date (Revision 1): FEB/25/2010